

**DIMAR**

Medical Devices

# *DimAir* Line

Intensive Care




***NIMV HELMET***

Single patient use Helmet for Non Invasive Mechanical Ventilation; it can be used also in CPAP and BI-PAP for the weaning from the ventilator, both in adults and children patients. It has 6 ISO connectors placed directly on the rigid ring, to offer more **EFFICACY, PRACTICALITY** and **COMFORT**. The version provided with pneumatically sealed zip, allows a quick and easy access to the patient.

- EFFICACY** The connection to the rigid ring reduces the delay of the trigger reaction time.
- PRACTICALITY** The double Right and Left symmetrical connectors allow a practical and easy connection with the equipment placed on any side.
- COMFORT** The connection between the Helmet and the tubes directly on the rigid ring makes the gas diffusion more comfortable and it reduces the inner noise.

## TECHNICAL FEATURES

- Opening zip with pneumatic seal\*
- Single patient use
- Made of biocompatible plastic material
- Soft and flexible straps with protective pads made of antibacterial CoolMax® material to reduce ulcerations
- 6 Right and Left symmetrical connectors
- 2 fixing eyelets on the top of the helmet for an additional/alternative connection
- 22M, 22F ISO connectors
- Latex Free



The “eyelet” additional connection prevents from ulceration under the ampits and improves inspiratory trigger reaction and the synchrony between patient and ventilator

## ADVANTAGES

- Effective, Practical and Comfortable
- Trigger reaction delay < 0,1”
- Quick access to the patient’s face
- Possibility to do nursing, feeding and weaning operations without removing the helmet
- Gases inlet and outlet can be inverted or connected on the same side; in this way the patients can also be placed on their side

## REFERENCES

CODE Standard	CODE * ZIP	SIZE	Neck circumference cm.	Internal volume** Lt. (except head)
600/5585	600/5586	XXS	23 ÷ 27	9
600/5228	600/5234	XS	27 ÷ 31	9
600/5229	600/5235	S	31 ÷ 35	9
600/5230	600/5236	M	35 ÷ 39	10
600/5231	600/5237	L	39 ÷ 43	10
600/5232	600/5238	XL	43 ÷ 47	11
600/5233	600/5239	XXL	47 ÷ 51	11
600/5580	600/5582	3XL	≥ 51	11



\*\*+/- dimensional tollerances, excluding the collar and the patient’s head

- Single pack
- Multiple pack with 5 pieces

## DIMAR

Via G. Galilei, 6 41036 Medolla -MO- I  
Tel: +39 0535 611336 Fax: +39 0535 611328  
Mail: dimar.srl@tiscali.it

