

A photograph of an elderly man with glasses and a flat cap, smiling broadly and holding a black cane. He is wearing a brown sweater with a dark argyle pattern. The background is a blurred outdoor setting. A large, semi-transparent version of his face is overlaid in the upper left corner.

**Mobility Aids**

# AIDS FOR THE ELDERLY AND DAILY LIVING

---

**TECHNOLOGY THAT CREATES  
SMILES**

### Palm Grip Cane (Right)



Product Code - PG003R

- Aluminium shaft, black finish height adjustable
- Ergonomic plastic moulded handle with rubber ferrule

### Palm Grip Cane (Left)



Product Code - PG003L

- Aluminium shaft, anodized finish, height adjustable
- Ergonomic plastic moulded grey handle with rubber ferrule

### Quad Cane (Small Base)



Product Code - DP001S

- Chrome Plated legs with aluminium shaft
- Large base for additional stability
- Height adjustable from 79 cm to 104 cm (31"-41")
- Rubber ferrules on all four legs
- This stick has four feet at the end of the cane shaft which provides a wider base offering more stability for the user

### Quad Cane (Large Base)



Product Code - DP001L

- Chrome Plated legs with aluminium shaft
- Large base for additional stability
- Height adjustable from 79 cm to 104 cm (31 inches - 41 inches)
- Rubber ferrules on all four legs
- This stick has four feet at the end of the cane shaft which provides a wider base offering more stability for the user

### Tripod Cane



**Product Code - DP002S**

- Chrome plated three legs with aluminium shaft
- Height adjustable from 79 cm to 104 cm (31" - 41")
- Rubber ferrules on all three legs

### Adjustable Cane



**Product Code - DP006**

- Aluminium shaft, height adjustable
- Anodized / powder coated
- Standard plastic moulded handle with rubber ferrule

### Double Adjustable Elbow Crutch with Standard Handle



**Product Code - EC001**

- Aluminium construction.
- PVC handgrips
- Height adjustment on leg and cuff
- Push-button adjustment makes it easy to adjust crutch height in 1" increments
- Grey rubber ferrules

### Double Adjustable Elbow Crutch with Soft Feel Ergo Handle



**Product Code - EC002**

- Crutches with anatomic comfort hand grips
- Grey rubber ferrules
- Push-button adjustment makes it easy to adjust crutch height in 1" increments

### Folding Walking Stick with Ergonomic Handle

**Product Code - FWS011**


- Folds in to 4 convenient parts for easy handling, black finish
- Attractive ergonomic handle combined with style and function
- Locking ring ensures added safety after adjustment
- Handle height adjusts from 33" - 37"
- Comes standard with wrist strap, lock silencer, carry bag
- Max user weight 120 kg



### Folding Walking Stick with Standard handle

**Product Code - FWS012**

- Folds in to 4 convenient parts for easy handling, black finish
- Attractive ergonomic handle combined with style and function
- Locking ring ensures added safety after adjustment
- Handle height adjusts from 33" - 37"
- Comes standard with wrist strap, lock silencer, carry bag
- Max user weight 90 kg



**Stainless Steel Grab Rail (12")**



**Product Code - GRSS12**

- Steel Tubing.
- Flanged with screw cover.

**Stainless Steel Grab Rail (18")**




**Product Code - GRSS18**

- Steel Tubing.
- Flanged with screw cover.

**Stainless Steel Grab Rail (24")**

**Product Code - GRSS24**

- Stainless steel 304 grade
- Anti-rust/ Anti corrosive
- Can be placed horizontally, vertically or diagonally
- Secure installation, Available Size 12" inches, 18" inches, & 24" inches



**Plastic Fluted Grab Rail (12")**



**Product Code - GRP12**

- Extruded rail
- Supplied with fixings
- Available in 12", 18", 21", 24" sizes
- Injection moulded end brackets with screw cover

**Plastic Fluted Grab Rail (18")**



**Product Code - GRP18**

- Extruded rail
- Supplied with fixings
- Available in 12", 18", 21", 24" sizes
- Injection moulded end brackets with screw cover

**Plastic Fluted Grab Rail (21")**



**Product Code - GRP21**

- Extruded rail
- Supplied with fixings
- Available in 12", 18", 21", 24" sizes
- Injection moulded end brackets with screw cover

**Plastic Fluted Grab Rail (24")**



**Product Code - GRP24**

- Extruded rail
- Supplied with fixings
- Available in 12", 18", 21", 24" sizes
- Injection moulded end brackets with screw cover

### Cane Seat (Foldable)

Product Code - CS007

- Convenient and light weight
- Strong aluminium construction with comfortable plastic molded seat & handle
- Tripod design with rubber ferrules



### Mobile Commode

Product Code - MC03

- Mobile commode can be used for toileting and bathing
- Mobile on 4 wheel but not for outside use
- PU seat, back, over seat & arm rest for comfort
- Detachable foot rest
- Fold - up arm rest







1701 North St Suite 40-2, Endicott,  
NY 13760, USA



[info@sanrai.com](mailto:info@sanrai.com)