





Product catalogue




SOMNOtouch™ RESP




Valid from 15th of June 2017





SOMNOtouch™ RESP: Sets








| Item no.: | Picture: | Description: |
|-----------|---|--|
| TOR100 |  | <p>SOMNOtouch™ RESP Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device, Belts(Abdomen/Thorax L), Effort sensor, SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), Bag, User manual, Software DOMINO light (2 licences)</p> |
| TOR100A |  | <p>Add-on device: SOMNOtouch™ RESP Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device, Belts (Abdomen/Thorax L), SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Bag</p> |
| TOR110 |  | <p>SOMNOtouch™ RESP (RIP) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device, 2 Inductive effort belts Abdomen/Thorax L + Body mount, inductive effort sensor (RIP), SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), Bag, User manual, Software DOMINO light (2 licences)</p> |
| TOR110A |  | <p>Add-on device: SOMNOtouch™ RESP (RIP) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size L + Body mount, Inductive effort sensor (RIP), SpO2 Finger clip Soft Silicone (Adult), 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Bag</p> |








| Item no.: | Picture: | Description: |
|-----------|---|--|
| TOR130 |  | <p>SOMNOtouch RESP (RIP) with thermal airflow recording Respiratory Screening</p> <p>Max. Recording time: 17 h</p> <p>12 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, Flow Thermistor, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size L + Body mount, Inductive effort sensor 2 channels (Thorax/Abdomen), Nasal/Oral flow sensor thermistor adult, SpO2 Finger clip Soft Silicone (adult), Retaining clip for simultaneous application of Thermistor and Nasal cannula, 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), Bag, User manual, Software DOMINO light (2 licences)</p> |
| TOR130A |  | <p>Add-on device: SOMNOtouch RESP (RIP) with thermal airflow recording Respiratory Screening</p> <p>Max. Recording time: 17 h</p> <p>12 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, Flow Thermistor, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size L + Body mount, Inductive effort sensor 2 channels (Thorax/Abdomen), Nasal/Oral flow sensor thermistor adult, SpO2 Finger clip Soft Silicone (adult), Retaining clip for simultaneous application of Thermistor and Nasal cannula, 1 pack Nasal cannulas (10 pcs.), CPAP adapter, Bag</p> |
| TOR105 |  | <p>SOMNOtouch™ RESP (pediatric) Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, Belt set Abdomen/Thorax size XS, Effort sensor, SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas ped. (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), Bag, manual, Software DOMINO light (2 licences)</p> |








| Item no.: | Picture: | Description: |
|-----------------------|---|---|
| <p>TOR105A</p> |  | <p>Add-on device : SOMNOtouch™ RESP (pediatric) Respiratory Screening</p> <p>Max. Recording time: 19 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse rate, Plethysmography, 2x Effort, Flow+Snoring (Nasal pressure), CPAP Pressure, Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, Belt set Abdomen/Thorax size XS, Effort sensor, SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas ped. (10 pcs.), CPAP adapter, Bag,</p> |
| <p>TOR115</p> |  | <p>SOMNOtouch RESP (RIP, pediatric) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Ind. effort belts Abdomen/Thorax size S + Body mount, Ind. effort sensor (RIP), SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas ped. (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), Bag, manual, Software DOMINO light (2 licences)</p> |
| <p>TOR115A</p> |  | <p>Add-on device: SOMNOtouch RESP (RIP, pediatric) Respiratory Screening</p> <p>Max. Recording time: 18 h</p> <p>11 channels: Body position, Movement, SpO2, Pulse, Pleth, 2x RIP Effort, Flow, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size S + Body mount, Inductive effort sensor (RIP), SpO2 Finger clip Soft Silicone (pediatric), 1 pack Nasal cannulas pediatric (10 pcs.), CPAP adapter, Bag</p> |








| Item no.: | Picture: | Description: |
|-----------|---|--|
| TOR135 |  | <p>SOMNOtouch RESP (RIP) with thermal airflow recording (pediatric) Respiratory Screening</p> <p>Max. Recording time: 17 h</p> <p>12 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, Flow Thermistor, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size S + Body mount, Inductive effort sensor 2 channels (Thorax/Abdomen), Nasal/Oral flow sensor thermistor pediatric, SpO2 Finger clip Soft Silicone (pediatric), Retaining clip for simultaneous application of thermistor and nasal cannula pediatric, 1 pack Nasal cannulas pediatric (10 pcs.), CPAP adapter, Dockingstation (incl. power supply and USB cable), Bag, User manual, Software DOMINO light (2 licences)</p> |
| TOR135A |  | <p>Add-on device: SOMNOtouch RESP (RIP) with thermal airflow recording (pediatric) Respiratory Screening</p> <p>Max. Recording time: 17 h</p> <p>12 channels: Body position, Movement, SpO2, Pulse, Plethysmography, 2x RIP Effort, Flow, Flow Thermistor, CPAP Pressure, Snoring (Nasal pressure), Patient marker</p> <p>Items included: Basic device SOMNOtouch RESP, 2 Inductive effort belts Abdomen/Thorax size S + Body mount, Inductive effort sensor 2 channels (Thorax/Abdomen), Nasal/Oral flow sensor thermistor pediatric, SpO2 Finger clip Soft Silicone (pediatric), Retaining clip for simultaneous application of thermistor and nasal cannula pediatric, 1 pack Nasal cannulas pediatric (10 pcs.), CPAP adapter, Bag</p> |

SOMNOtouch™ RESP: Options for extended analysis










| Item no.: | Picture: | Description: |
|-----------|---|---|
| TOS010S |  | 3-channel ECG option for SOMNOtouch NIBP and RESP Items included: 3 channel ECG sensor, 1 pack of disposable electrodes, short instruction manual |
| TOS011S |  | 1-channel ECG option for SOMNOtouch RESP Items included: 1 channel ECG sensor, 1 pack of disposable electrodes, short instruction manual |
| TOS012S |  | 1-channel ECG option pediatric for SOMNOtouch RESP Items included: 1 channel ECG sensor pediatric, 1 pack of disposable electrodes, short instruction manual |
| TOS014S |  | 3-channel ECG option pediatric for SOMNOtouch RESP and NIBP Items included: 3-channel ECG sensor pediatric, 1 pack of disposable electrodes, short instruction manual |
| TOS016S | <i>Available soon</i> | 1-channel ECG/IC-EMG option for SOMNOtouch RESP Items included: 1-channel ECG/IC-EMG sensor, 2 pack of disposable electrodes, short instruction manual +++ Available soon +++ |
| TOS041 |  | 1-channel EEG option for SOMNOtouch RESP Items included: 1-channel EEG sensor, Activation of EEG option, 1 pack of disposable electrodes |
| TOS051 |  | Simplified Sleep Stage Evaluation (ambulatory) for SOMNOtouch RESP Items included: EEG / EOG Combi sensor (3 channels) with integrated light sensor, Activation of EEG option, 1 pack of disposable electrodes |
| TOS095 |  | AASM Headbox for SOMNOtouch RESP Items included: AASM Headbox (4 EEG, 2 EOG, 1 EMG/ECG, Impedance check), Activation of AASM Sleep Option, 10 GRASS gold plated EEG cup electrodes, EC2 Electrode cream (self-adhesive), NUPREP Abrasive skin prepping gel |






| Item no.: | Picture: | Description: |
|-----------|---|---|
| TOS097 |  | <p>R&K Headbox for SOMNOtouch RESP</p> <p>Items included: R&K Headbox (4 EEG, 2 EOG, 1 EMG/ECG, Impedance check), Activation of R&K Sleep Option, 10 GRASS gold plated EEG cup electrodes, EC2 Electrode cream (self-adhesive), NUPREP Abrasive skin prepping gel</p> |
| TOS091 |  | <p>6-channel EEG option for SOMNOtouch RESP</p> <p>Items included: 6-channel EEG module (6 EEG, 1 EMG/ECG, Impedance control), Activation of EEG option, 10 GRASS gold plated EEG cup electrodes (10 mm with hole, Cable length: 76 cm), EC2 Electrode cream (self-adhesive), NUPREP Abrasive skin prepping gel</p> |
| TOS081 |  | <p>Set for inductive effort measurement (RIP) for SOMNOtouch RESP</p> <p>Items included: Inductive effort sensor 2 channels (Thorax/Abdomen), Inductive effort belts (2 pieces, size L), Body Mount SOMNOtouch for inductive effort belts</p> |
| TOS110S |  | <p>2-channel PLM Option for SOMNOtouch RESP</p> <p>Items included: 2-channel PLM sensor, 1 pack of disposable electrodes, short instruction manual</p> |
| TOS140S |  | <p>Effort + 1PLM Combi Option for SOMNOtouch RESP</p> <p>Items included: Effort + 1PLM Combi Sensor , 1 pack of disposable electrodes, short instruction manual</p> |
| TOS170S |  | <p>Effort + 2PLM Combi Option for SOMNOtouch RESP</p> <p>Items included: Effort + 2PLM Combi Sensor, 1 pack of disposable electrodes, short instruction manual</p> |
| TOS175S |  | <p>Effort + 2PLM + 1 channel ECG Combi Option for SOMNOtouch RESP</p> <p>Items included: Effort + 2 PLM + 1-channel ECG Combi Sensor, 2 pack of disposable electrodes, short instruction manual</p> |

| Item no.: | Picture: | Description: |
|--------------|---|--|
| TOS191 |  | <p>Set for inductive effort measurement (RIP) and thermal airflow recording for SOMNOtouch RESP</p> <p>Items included: Inductive effort sensor 2 channels (Thorax/Abdomen), Nasal/Oral flow sensor thermistor adult,- Retaining clip for simultaneous application of Thermistor and Nasal cannula, Inductive effort belts (2 pieces, size L), Body Mount SOMNOtouch for inductive effort belts</p> |
| TOS192 |  | <p>Set for inductive effort measurement (RIP) and thermal airflow recording for SOMNOtouch RESP (pediatric)</p> <p>Items included: Inductive effort sensor 2 channels (Thorax/Abdomen), Nasal/Oral flow sensor thermistor pediatric,- Retaining clip for simultaneous application of Thermistor and Nasal cannula pediatric, Inductive effort belts (2 pieces, size S), Body Mount SOMNOtouch for inductive effort belts</p> |
| TOS151 |  | <p>Pneumotachograph Set for SOMOtouch RESP</p> <p>Items included: tube for the connection of CPAP/BiPAP devices, analyses the quantitative breathing and leakage flow</p> |
| TOS201 |  | <p>Pneumotachograph Set for SOMOtouch RESP</p> <p>incl. tube for the connection of CPAP/BiPAP devices, analyses the quantitative breathing and leakage flow</p> <p>+++ Available soon as successor of TOS151 +++</p> |
| VID98A |  | <p>LAN Infrared Video Camera Home Sleep</p> <p>H.264 Video Stream, Integrated IR illumination, Microphone, Resolution 800x600 pixels, Power supply, USB stick 16 GB for data storage, USB to Ethernet Adapter, Lan cable, tripod, Carrying bag, Software for synchronized recording and archiving</p> |
| TOS301 |  | <p>Blood Pressure Option for SOMNOtouch RESP</p> <p>Items included: 1-channel ECG sensor for SOMNOtouch RESP, Activation of BP option (based on serial number), 1 Pack of Adhesive Electrodes</p> |
| KOMP112 D |  | <p>Option Signal Check via Tablet PC (Bluetooth) including SOMNOmedics Android Application</p> <p>to ensure the functionality of the live data transmission, a SOMNO HD, SOMNOscreen BT Series or a SOMNOtouch is needed</p> |














| Item no.: | Picture: | Description: |
|----------------|---|---|
| KOMP111 |  | <p>Option Signal Check via smartphone including SOMNOmedics Android application</p> <p>Signal check via e-mail for ambulatory polygraphy and polysomnography</p> <p>+++ This option requires a valid smartphone contract +++</p> |
| SSC088 |  | <p>Licence Server</p> <p>Enables filing of the software licences at a central server. Because of the central administration an easy distribution of the software to the single working stations is possible by software packet assembling. This software packing will contribute to easy software maintenance. Additionally, users can flexibly log on to different working stations, which facilitates scoring of recordings at individual working stations as well as an easy embedding of external scorers.</p> |
| SSC091 |  | <p>DOMINO light Software option: extended Data Base for measurement results</p> |
| SW005 |  | <p>DOMINO light Software option: Online signal display and recording for SOMNOtouch RESP</p> |
| SW006 |  | <p>DOMINO light Software option: Online signal display and recording (incl. Bluetooth-Stick) for SOMNOtouch RESP</p> |
| SW010 |  | <p>Schiller Darwin OEM, ECG-Arrhythmia Analysis for DOMINO (light)- valid for one PC</p> |
| SOW105 |  | <p>DOMINO light Additional Software License</p> |

SOMNOtouch™ RESP: Sensors/Accessories

| Item no.: | Picture: | Description: |
|-----------|---|---|
| TOS020 |  | Effort Sensor für SOMNOtouch RESP |
| TOS060 |  | Soft Silicone Finger Sensor for SOMNOtouch RESP 1.3 m, for adults |
| TOS062 |  | Soft Silicone Finger Sensor for SOMNOtouch RESP 1.0 m, pediatric |
| TOS064 |  | Soft Silicone Finger Sensor for SOMNOtouch RESP (size: medium) 1.3 m, for adults |
| TOS070 |  | SoftWrap Finger Sensor for SOMNOtouch RESP 1.3 m, for adults and pediatric use |
| TOS080 |  | Inductive effort sensor (RIP) for SOMNOtouch RESP 2 channels (Thorax/Abdomen) |
| TOS085 |  | Inductive effort sensor (RIP) pediatric for SOMNOtouch RESP, 2 channels (Thorax/Abdomen) 1m cable |
| TOS120 |  | Analogue Optocoupler Connecting external devices to the SOMNOtouch Requires technical data sheet of the external devices! (An additional charge for the adaption depending on the needed equipment will be added) |
| TOS121 |  | 4-channel analogue Optocoupler Connecting external devices to the SOMNOtouch Requires technical data sheet of the external devices! (An additional charge for the adaption depending on the needed equipment will be added) |




| Item no.: | Picture: | Description: | |
|-----------|---|---|-------------|
| TOS190 |  | RIP-Thermistor Combi Sensor for SOMNOtouch RESP | |
| NGS151 |  | Nasal/Oral flow sensor thermistor adult | |
| NGS152 |  | Nasal/Oral flow sensor thermistor pediatric | |
| TOS210 |  | Bruxism sensor for SOMNOtouch RESP +++ Available soon +++ | |
| TOS902 |  | Docking station SOMNOtouch incl. power supply and USB cable | for |
| TOS903 |  | 3-port Docking station for SOMNOtouch supply and USB cable | incl. power |
| SOW080 |  | Bag for SOMNOWatch plus and SOMNOtouch | |

SOMNOtouch™ RESP: Belts






| Item no.: | Picture: | Description: |
|-----------|---|--|
| TOS799 |  | Effort belt for effort sensor SOMNOtouch RESP size XS, circumference: ca. 50 - 75 cm, fleece with quick release and velcro fastener |
| TOS800 |  | Effort belt for effort sensor SOMNOtouch RESP size S, circumference: ca. 65 - 88 cm, fleece with quick release and velcro fastener |
| TOS801 |  | Effort belt for effort sensor SOMNOtouch RESP size M, circumference: ca. 74 - 109 cm, fleece with quick release and velcro fastener |
| TOS802 |  | Effort belt for effort sensor SOMNOtouch RESP size L, circumference: ca. 86 - 134 cm, fleece with quick release and velcro fastener |
| TOS803 |  | Effort belt for effort sensor SOMNOtouch RESP size XXL, circumference: ca. 99 - 159 cm, fleece with quick release and velcro fastener |
| TOS804 |  | Effort belt for SOMNOtouch RESP, size XS, Circumference: ca. 50 - 75 cm, fleece with quick release and velcro fastener |
| TOS805 |  | Effort belt for SOMNOtouch RESP size S, Circumference: ca. 65 - 88 cm, fleece with quick release and velcro fastener |
| TOS806 |  | Effort belt for SOMNOtouch RESP size M, Circumference: ca. 74 - 109 cm fleece with quick release and velcro fastener |
| TOS807 |  | Effort belt for SOMNOtouch RESP size L, Circumference: ca. 86 - 134 cm, fleece with quick release and velcro fastener |
| TOS808 |  | Effort belt for SOMNOtouch RESP size XXL, Circumference: ca. 99 - 159 cm, fleece with quick release and velcro fastener |
| TOS810 |  | Body Mount SOMNOtouch RESP for inductive effort belts, Width 38 mm, Length 15 cm |
| TOS811 |  | Application Mount SOMNOtouch RESP |
| TOS812 |  | Body Mount SOMNOtouch RESP for inductive disposable effort belts, Width 38 mm, Length 18 cm |

| Item no.: | Picture: | Description: |
|-----------|---|--|
| TOS813 |  | Application Mount with retaining clip for SOMNOtouch RESP |
| SEN660 |  | Inductive Effort Belts - Size S Pediatric, (2 pcs./Pack), Circumference: ca. 50 - 75 cm |
| SEN661 |  | Inductive Effort Belts - Size M Pediatric and Adult, (2 pcs./Pack), Circumference: ca. 60 - 90 cm |
| SEN662 |  | Inductive Effort Belts - Size L Adult, (2 pcs./Pack), Circumference: ca. 80 - 120 cm |
| SEN663 |  | Inductive Effort Belts - Size XXL Adult, (2 pcs./Pack), Circumference: ca. 100 - 150 cm |
| SEN664 |  | Inductive Effort Belts - Size XS Pediatric, (2 pcs./pack) Circumference: ca. 35-50 cm |
| SEN680 |  | Disposable inductive Effort Belt - Size S Circumference: ca. 40-70 cm, (4 pcs./ Pack) ++ May not be available in all countries due to regulatory restrictions +++ |
| SEN681 |  | Disposable inductive Effort Belt - Size M Circumference: ca. 60-90 cm, (4 pcs./ Pack) ++ May not be available in all countries due to regulatory restrictions +++ |
| SEN682 |  | Disposable inductive Effort Belt - Size L Circumference: ca. 90-120 cm, (4 pcs./ Pack) ++ May not be available in all countries due to regulatory restrictions +++ |
| SEN683 |  | Disposable inductive Effort Belt - Size XXL Circumference: ca. 120-150 cm, (4 pcs./ Pack) ++ May not be available in all countries due to regulatory restrictions +++ |

SOMNOtouch™ RESP: Creams

| Item no.: | Picture: | Description: |
|-----------|---|---|
| SEN030 |  | Nuprep Abrasive skin prepping gel, 114 g per tube, 3 pcs./Pack |
| SEN032 |  | Ten20 EEG Conductive Adhesive Paste, 114 g tube, 3 pcs./Pack |
| SEN034 |  | EC2 Adhesive electrode cream, water-based, 100 g per tube, 10 pcs./Pack |






SOMNOtouch™ RESP: EEG-Electrodes

| Item no.: | Picture: | Description: |
|-----------|---|---|
| SEN018 |  | Connecting cables for disposable electrodes Snap fastener electrode cables (∅ 1.5 mm Safety pins), Cable length 80 cm, assorted colours, 10 pcs./Pack |
| F-E5GH-24 |  | GRASS Gold Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 61 cm, 10 pcs./Box |
| F-E5GH-30 |  | GRASS Gold Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 76 cm, 10 pcs./Box |
| F-E5GH-48 |  | GRASS Gold Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 122 cm, 10 pcs./Box |
| F-E5SH-48 |  | GRASS Silver Cup Electrodes 10 mm with hole 2 mm, assorted colours, Cable length 122 cm, 10 pcs./Box |

SOMNOtouch™ RESP: Electrodes and consumables

| Item no.: | Picture: | Description: |
|-----------|---|--|
| SEN002 |  | Single use electrode N (recommended for ECG) Single use electrode with decentral snap fastener, 25 pcs./Pack |
| SEN004 |  | Adhesive Electrode (recommended for ECG) Skin-friendly Santara-microporous fleece textile, 30 pcs./Pack |
| SEN006 |  | Blue Sensor Disposable Electrodes (recommended for ECG) Snap fastener, Diameter 6.8 cm, 25 pcs./Pack |
| SEN010C |  | Disposable Electrodes (recommended for EOG, EEG) Silver-Silver chloride, self-adhesive, Snap fastener, Diameter 2.5 cm, 30 pcs./Pack |
| SEN011 |  | Disposable Electrodes Silver-Silver chloride, self adhesive, Snap fastener, Diameter 3 cm, 30 pcs./Pack |
| SEN012 |  | Disposable Electrodes (recommended for PLM and ECG) Silver-Silver chloride, self-adhesive, Snap fastener, 2.8 x 4.5 cm, 30 pcs./Pack |
| SEN012A |  | Disposable Electrodes (recommended for ECG and PLM) Silver-Silver chloride, self adhesive, Snap fastener, 2.8 x 4.5 cm, 30 pcs./Pack ++ May not be available in all countries due to regulatory restrictions +++ |
| SEN013A |  | Disposable Electrodes (recommended for Intercostal EMG and PLM) 20 pcs./ Pack ++ May not be available in all countries due to regulatory restrictions +++ |

| Item no.: | Picture: | Description: |
|-----------|---|---|
| SEN020 |  | Disposable Electrode with connecting cable Clear gel and Safety plug (colour coded) 1.5 mm, Cable length 60 cm, Diameter 2.2 cm, 15 pcs./Pack |
| SEN020A |  | Disposable Electrode with connecting cable Clear gel and Safety plug (colour coded) 1.5 mm, Cable length 120 cm, Diameter 2.2 cm, 15 pcs./Pack |
| SEN020B |  | Disposable Electrode with connecting cable Clear gel and Safety plug (colour coded) 1.5 mm, Cable length 160 cm, Diameter 2.2 cm, 15 pcs./Pack |
| TOS075 |  | Wrap tapes for SoftWrap Fingerclip 12 pcs./Box, Single to use, highly flexible and ultra-soft wrap tape |
| SEN115 |  | Nasal Cannula HC (High Comfort) Adult, Nasal Pressure transducer with LuerLock connector male, Tube length 60 cm, 10 pcs./Pack |
| SEN115A |  | Nasal Cannula Adult, Nasal Pressure transducer with LuerLock connector male, Tube length 60 cm, 100 pcs./Pack ++ May not be available in all countries due to regulatory restrictions ++ |
| SEN116B |  | Nasal Cannula Pediatric Nasal Pressure transducer with LuerLock connector male, Tube length 40 cm, 10 pcs./Pack |
| SEN117 |  | Nasal Cannula Neonatal, Nasal Pressure transducer with LuerLock connector male, Tube length 210 cm, 25 pcs./Pack |
| SEN118 |  | Nasal/Oral Cannula for Pressure Sensor Adult, Nasal and Oral Pressure transducer with LuerLock connector male, Tube length 60 cm, 25 pcs./Pack |

| Item no.: | Picture: | Description: |
|-----------|---|---|
| SEN118A |  | <p>Nasal/Oral Cannula for Pressure Sensor Adult, Nasal and Oral Pressure transducer with LuerLock connector male, Tube length 60 cm, 100 pcs./Pack</p> <p>++ May not be available in all countries due to regulatory restrictions +++</p> |
| SEN121 |  | <p>Nasal Cannula (divided) Adult, for Airflow Pressure/Snore monitoring and simultaneous gas sampling, tube length 210 cm, 10 pcs./pack</p> |
| SEN903 |  | <p>Nasal Cannula (CO2) Adult, for Capnometry, Tube length 210 cm, 10 pcs./Pack</p> |
| SEN921 |  | <p>Silicone Tube for CPAP connection with LuerLock</p> |
| SEN924 |  | <p>CPAP adapter</p> |

****** PCs/Notebooks are available in different configurations on request*****