



# Homecare Product Catalog

# CONTENT

## CPAP

	Page
Masks and Replacement Parts	3-5
Tubing	6-7
Filters	8-10
Respironics	8
ResMed	9
DeVilbiss	10
Resvent	10
Humidifier Chambers	10
Chin straps	11-12
Headgear	12
Cash Items	13-20
Cleaners and Accessories	13-14
Batteries	15-16
Mask Liners and Soft Headgear Covers	17
Self Care	18
Connectors, Adapters, and Inline Filter	18
Pillows	19
Storage	19
Pressure Meter	20
Tubing Wraps	20

## Oxygen

	Page
Pulse Oximeter	21
Firesafe	22
Masks	23
Cannula	24-25
Tubing	25
Cylinder Carriers	26-27
Connectors and Adapters	28
Accessories	29-30
Cleaners	30
Bacteria and HEPA Filters	31
Concentrator Cabinet Filters	32-33
Caire/AirSep/SeQual	32
DeVilbiss/Drive	32
Respironics	33
Nidek	33
Jumao	33

## Other

	Page
Aerosol	35
Handheld Nebulizers	36
Nebulizer Kits	37
Nebulizer and Anesthesia Filters	38
Suction	39-40
Closed Suction	41-42
Tracheostomy	43
Ventilator Filters, Masks, and Connectors	44
Ventilator Circuits	45-46
StrongArm Comfort Cane	47
Anesthesia	48-50
Anestand	48-49
GuideBlade	50

# SUNSET HEALTHCARE SOLUTIONS

## Sol Full Face CPAP Mask

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM205S</b>
Medium	<b>CM205M</b>
Large	<b>CM205L</b>

**NEW**

Ask us about our  
30-day mask guarantee

### Replacement Cushions

Size	Item#
Small	<b>CUCM205S</b>
Medium	<b>CUCM205M</b>
Large	<b>CUCM205L</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031



### A premium mask without the premium price

#### Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

#### Quiet vent


Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

#### Quick release connector

Rapidly and easily separates or connects mask and the tube

#### Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment



Sol comes in attractive retail packaging



## Ray Pediatric CPAP/NIV Masks

### Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	<b>CM061S</b>
Medium	<b>CM061M</b>
Large	<b>CM061L</b>
Extra Large	<b>CM061XL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM061S</b>
Medium	<b>CUCM061M</b>
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031



**NEW**

### Nasal Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	<b>CM060-KITS</b>
L&XL Cushions	<b>CM060-KITL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM060S</b>
Medium	<b>CUCM060M</b>
Large	<b>CUCM060L</b>
Extra Large	<b>CUCM060XL</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032



• Comparable to Nonny masks



# SUNSET HEALTHCARE SOLUTIONS

## Sunset Deluxe Nasal Pillow CPAP Mask

Item#

**CM036**

Mask with Headgear and Small, Medium, and Large Pillows

**HG036**

Headgear Only

HCPCS#

MASK - A7034

HEADGEAR - A7035

PILLOWS - A7033



### Replacement Pillows

Size	Item#
Small	<b>CM036SP</b>
Medium	<b>CM036MP</b>
Large	<b>CM036LP</b>



### Mask without Headgear

**CM036NH**

- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

## Sunset Clearsight Deluxe Nasal Mask

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM110S</b>
Medium	<b>CM110M</b>
Large	<b>CM110L</b>

**HG003** Headgear - Fits Most

**HG003XL** Headgear - XL

### Replacement Cushions

Size	Item#
Small	<b>CUCM110S</b>
Medium	<b>CUCM110M</b>
Large	<b>CUCM110L</b>

HCPCS#

MASK - A7034

HEADGEAR - A7035



- Clearsight design provides unobstructed view

Mask Series	Mask Type	Replaceable Cushion	Forehead Support	Flex Tube Included	Headgear Included
<b>CM005</b>	Full Face	✓		✓	✓
<b>CM105</b>	Full Face	✓	✓ (Adjustable)		✓
<b>CM006</b>	Nasal	✓	✓	✓	✓
<b>CM106</b>	Nasal	✓	✓ (Adjustable)		✓
<b>CM007</b>	Full Face				✓
<b>CM008</b>	Nasal				✓
<b>CM110</b>	Nasal	✓			✓
<b>CM036</b>	Nasal Pillow	✓		✓	✓

## Sunset Classic Full Face CPAP Mask

### Mask with Headgear

Size	Item#
Small	<b>CM007S</b>
Medium	<b>CM007M</b>
Large	<b>CM007L</b>
Extra Large	<b>CM007XL</b>

**HG003** Headgear - Fits Most

**HG003XL** Headgear - XL

**CMSIZE02**

Classic Sizing Guide



HCPCS#  
MASK - A7030  
HEADGEAR - A7035



**1/2 THE WEIGHT  
OF LEADING MASKS**

- Conventional forehead pad removed
- to create a clear line of sight
- Minimalist design and ultra-lightweight silicone
- Quick-release clips make it easy to remove
- 360° dual-swivel elbow port allows for free movement
- One-piece mask and cushion design is durable and easy

## Sunset Classic Nasal CPAP Mask

### Mask with Headgear

Size	Item#
Medium	<b>CM008M</b>
Large	<b>CM008L</b>

**HG003** Headgear - Fits Most

**HG003XL** Headgear - XL

**CMSIZE02**

Classic Sizing Guide



HCPCS#  
MASK - A7034  
HEADGEAR - A7035

- Great value for cash patients or those with lower insurance reimbursement
- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

# SUNSET HEALTHCARE SOLUTIONS

## Sunset Deluxe Full Face CPAP Masks

### Mask with Flex Tube, Headgear, and Replaceable Cushion

Size	Item#
Small	<b>CM005S</b>
Medium	<b>CM005M</b>
Large	<b>CM005L</b>

Unobtrusive design offers a clear line of sight


**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPSC#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031

### Mask without Headgear

Size	Item#
Small	<b>CM005SNH</b>
Medium	<b>CM005MNH</b>
Large	<b>CM005LNH</b>



### Replacement Cushions

Size	Item#
Small	<b>CUCM005S</b>
Medium	<b>CUCM005M</b>
Large	<b>CUCM005L</b>

- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone replaceable cushion increases comfort while maintaining optimal seal

### Mask with Adjustable Forehead T, Headgear, and Replaceable Cushion

Size	Item#
Small	<b>CM105S</b>
Medium	<b>CM105M</b>
Large	<b>CM105L</b>



## Sunset Deluxe Nasal CPAP Masks

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM006S</b>
Medium	<b>CM006M</b>
Large	<b>CM006L</b>

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPSC#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032

### Mask without Headgear

Size	Item#
Small	<b>CM006SNH</b>
Medium	<b>CM006MNH</b>
Large	<b>CM006LNH</b>



### Replacement Cushions

Size	Item#
Small	<b>CUCM006S</b>
Medium	<b>CUCM006M</b>
Large	<b>CUCM006L</b>

### Mask with Headgear, Adjustable Forehead T, and Replaceable Cushion

Size	Item#
Small	<b>CM106S</b>
Medium	<b>CM106M</b>
Large	<b>CM106L</b>



# CPAP TUBING

## Durable CPAP Tubing — Gray



Item#

**TUB06**

6ft

19mm with 22mm cuffs

HCPCS# A7037

- Sunset's most popular item
- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free


### All Lengths

<b>TUB0018</b>	18in
<b>TUB002</b>	2ft
<b>TUB004</b>	4ft
<b>TUB06</b>	6ft
<b>TUB006.5</b>	6ft 6in — Compatible with ResMed tube wraps
<b>TUB007</b>	7ft
<b>TUB008</b>	8ft
<b>TUB009</b>	9ft
<b>TUB010</b>	10ft

Compatible with Sunset's new quick release connector!

Item#

**CAP1016**



## Slim CPAP Tube

Item#

**TUB006SS**

6ft x 15mm with  
22mm cuffs

HCPCS# A7037

### Highly durable

Strain relief cuff decreases separation and tearing

### Virtually no memory

Designed to be the most flexible tube available

### Minimal weight

One of the lightest tubes on the market at 2.6 ounces

### Great value


Affordable and sustainable



## Pressure Port Tubing

Item# **TUB006PP**

6ft with Pressure Sensor Lines  
22mm Cuffs • HCPCS# A7037



- Built-in sensor line for use with the Puritan Bennett: GoodKnight 420S, 420E and 425 and KnightStar 330

# HEATED CPAP TUBING

## Heated Tube for Respironics

Item#


**TUB06HPR**

Heated Tube for  
Philips Respironics  
Dreamstation Series  
and System One

OEM# HT15

HPCPS# A4604

- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off



**TUBHC900ICON208**

Heated Tube for Fisher & Paykel ICON

OEM# 900ICON208

## Hybernite Heated CPAP Tube 6ft

Item#

**TUB006ROC**

*MOQ may apply*

**NEW**

- Copper wire embedded in the tubing wall generates uniform heat and maintains the air temperature inside the tube to prevent rainout.
- Stand-alone kit can be used in the sleep lab or at home with any CPAP or BiPAP machine with humidifier.
- Includes heated circuit tube and independent power supply.

Compatible with machines that use standard CPAP Tubing.  
Not compatible with machines that use slim 15mm tubing  
such as newer Respironics and ResMed machines.



**TUB006HBT**

Hybernite Replacement Heated CPAP Tube 6ft



## DreamStation 2

**NEW**



Item# **CF1012**  
Disposable Filter  
**6/pack**

OEM# 1142831

HCPCS# A7038

**CF1012-1** 1/pack

**CF1012-2** 2/pack

**CF1012-3** 3/pack

**CF1012-50** 50/pack

Item# **CF1012F**  
Reusable Filter  
**2/pack**

OEM# 1142687

HCPCS# A7039

**CF1012F-1** 1/pack

## DreamStation



Item# **CF1010**  
Disposable Filter  
**6/pack**

OEM# 1122519

HCPCS# A7038

**CF1010-1** 1/pack

**CF1010-2** 2/pack

Item# **CF1010F**  
Reusable Filter  
**2/pack**

OEM# 112446 (1-pack)

HCPCS# A7039

**CF1010F-1** 1/pack

**CF1010F-6** 6/pack

## DreamStation Go Filters



Item# **CF101130D**  
30-Day Disposable Filter  
**6/pack**

OEM# 1133744

HCPCS# A7038

Item# **CF101190D**  
90-Day Reusable Filter  
**2/pack**

OEM# 1133743

HCPCS# A7039



## AirSense 11



Item#  
**CF2111**

Disposable Filter  
**3/pack**

HCPCS# A7038

**CF2111-1** 1/pack OEM# 39300

**CF2111-2** 2/pack

Item#

**CF2111HA-3**

Hypoallergenic Filter  
**3/pack**

1 1/4" x 1 3/16"

HCPCS# A7038

**CF2111HA-1** 1/pack

**CF2111HA-2** 2/pack

**NEW**

## S9 and AirSense 10 Series



Item#

**CF2107**

Disposable Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107-1** 1/pack OEM# 36850

**CF2107-2** 2/pack OEM# 36851

**CF2107-50** 50/pack OEM# 36853



Item#

**TUB06S9COMBO**

S9 Filter and CPAP Tube Kit

6 - CF2107

1 - 6ft CPAP Tube



**KIT**

Item#

**CF2107HA**

Hypoallergenic Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107HA-1** 1/pack OEM# 36855

**CF2107HA-2** 2/pack OEM# 36856

**CF2107HA-50** 50/pack

OEM# 36858



## AirMini



Item#

**CF2108**


Disposable Filter

**2/pack** • 1 1/4" x 1 3/16"

OEM# 38815

HCPCS# A7038

## S8 Series



Item#

**CF2105**

Disposable Filter

**3/pack** • 1 3/8" x 1 3/8"


OEM# 33916

HCPCS# A7038

**CF2105-1** 1/pack OEM# 33915

**CF2105-2** 2/pack OEM# 33973

**CF2105-50** 50/pack OEM# 33918



Item#

**CF2105HA**

Hypoallergenic Filter

**3/pack** • 1 3/8" x 1 3/8"

OEM# 33920

HCPCS# A7038

**CF2105HA-1** 1/pack OEM# 33919

**CF2105HA-2** 2/pack OEM# 33975

**CF2105HA-50** 50/pack OEM# 33922

# DeVILBISS

## IntelliPAP



Item# **CF5004**  
Ultrafine Filter  
**4/pack** • 1<sup>3</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
OEM# DV51D-603  
HCPCS# A7038

Item# **CF5004F**  
Foam Filter  
**4/pack** • 1<sup>5</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>16</sub>"  
OEM# DV51D-602  
HCPCS# A7039

## Horizon LT



Item# **CF5003**  
Ultrafine Filter  
**4/pack** • 4<sup>1</sup>/<sub>4</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
OEM# 8000D-603  
HCPCS# A7038

Item# **CF5003F**  
Foam Filter  
**4/pack** • 4<sup>1</sup>/<sub>4</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
OEM# 8000D-602  
HCPCS# A7039

# RESVENT

## iBreeze



HCPCS# A7038

Ultrafine Filter  
**CF9901-1** 1/pack  
**CF9901-2** 2/pack  
**CF9901-6** 6/pack  
**CF9901-50** 50/pack

# HUMIDIFIER WATER CHAMBERS

## Humidifier Chamber for ICON Series



Item# **HUMHC900ICON200**  
Replacement Chamber for the  
Fisher & Paykel ICON Series  
**Each** OEM# 900ICON200 • HCPCS# A7046

Compatible with:

- ICON Premo
- ICON Novo
- ICON Auto



## Humidifier Chamber for PR System One



Item# **HUM1136502**  
Humidifier Chamber for PR System  
One 60 Series  
**Each** HCPCS# A7046

## Humidifier Chamber for DreamStation



Item# **HUM1122520**  
Humidifier Chamber  
for DreamStation  
**Each** HCPCS# A7046

**NEW**

# SUNSET CHINSTRAPS

All chinstraps are latex-free • Hand wash only • Lay flat to dry

## Universal Neoprene Chinstrap



Item#

**CS001**

Universal Chinstrap

HCPCS# A7036



- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

## Basic Chinstrap



Item# **CS002**

Compare to  
Respironics Basic

OEM# 302175

HCPCS# A7036



- One size fits most
- Attaches to the top strap of headgear
- Cups the chin for comfort
- Works best with any nasal mask

## Adjustable Neoprene Chinstrap



Size

Item#

6in

**CS003S**

7in

**CS003R**

8in

**CS003M**

9in

**CS003L**

10in

**CS003XL**

Compare to Puresom Ad-  
justable Ruby

OEM# TMS-09 ADJ

HCPCS# A7036

- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head.
- Great for hard-to-fit patients

For fitting, measure along one side of the head, as indicated by the dashed line above.

## Deluxe Chinstrap



Size

Item#

26in

**CS004S**

28in

**CS004L**

Compare to  
Respironics Deluxe

OEM# 302425

HCPCS# A7036



- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

## Comfort Chinstrap



Item# **CS005**

Compare to  
ResMed Sullivan

OEM# 16015

HCPCS# A7036



- One size fits most
- Easy to use —Attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

## Premium Chinstrap



Item# **CS006**

Compare to  
Respironics Premium

OEM# 1012911

HCPCS# A7036




- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

# SUNSET CHINSTRAPS AND HEADGEAR

All chinstraps and headgear are latex-free • Hand wash only • Lay flat to dry

## Ruby Style Chinstrap



Size	Item#
21.5-22.5in	<b>CS007S</b>
22-23in	<b>CS007M</b>
23-24in	<b>CS007L</b>

Compare to Puresom Ruby  
OEM# TMS-09 • HCPCS# A7036

- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

## Halo Style Chinstrap



Item# **CS025**

Compare to BreatheWear  
Halo Chinstrap

OEM# H3407  
HCPCS# A7036

## Ruby Style Adjustable Chinstrap



Size	Item#
Regular	<b>CS007-ADJ</b>
Extra Large	<b>CS007-ADJXL</b>

HCPCS# A7036

- Fully adjustable
- Super soft materials for extra comfort



## Universal Mesh Headgear



Item# **HG001**

Compare to Respiroics  
SimpleStrap Headgear

OEM# 302328  
HCPCS# A7035

- One size fits most
- Can be used with most 4-point nasal or full face masks



## Universal Headgear



Size	Length	Item#
Regular	22.5in	<b>HG003</b>
Extra Large	24.5in	<b>HG003XL</b>

HCPCS# A7035

- Can be used with most 4-point nasal or full face masks



# CLEANERS AND ACCESSORIES

## Sunset Cleaning Wipes — Flowpack

Item#  
**CAP1003S**

**12/case**  
Flowpack of 64 wipes

**CAP1003S-80**

**12/case**  
Flowpack of 80 wipes



- Ingredients are safe and gentle for daily cleaning of all types of CPAP and oxygen masks and equipment

## CleanSmart CPAP Disinfectant Spray

Item#  
**CAP1018-L**  
16oz • **6/case**

**CAP1018-S**  
2oz • **9/case**

- Kills 99.9% of viruses & bacteria
- Included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2
- Active ingredient Hypochlorous Acid does not leave a chemical residue
- EPA Registration Number 89896-2
- Formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment



## Citrus II Cleaning Wipes

Item# **CAP1003**  
Canister of 62 wipes  
**12/case**



## Citrus II Mask Cleaning Spray

Item# **CAP1001**  
1.5oz • **24/case**

Item# **CAP1002**  
8oz • **12/case**

- All-natural cleanser for daily mask cleaning
- Available in 1.5oz travel size or 8oz non-aerosol spray cans



## Cleaning Wipes — Single Packs

**CAP1005**  
Citrus II Box of 12 single-pack wipes  
**12/case**

- Convenient individually wrapped wipes are great for patient samples or travel use




# CLEANERS AND ACCESSORIES

## CPAP Tube Cleaning Brush


Size	Item#
Standard 6ft	<b>CAP1013</b>
Slim 6ft	<b>CAP1013SS</b>

- Brushes all the way through a 6ft CPAP tube without leaving debris in the middle
- Stainless steel wire and soft nylon bristles



## Filter Replacement Kit for SoClean 2

Item# **CAP1007S-KIT**  
Each



- With filter, check valve, retail packaging and instruction manual
- Replace every 6 months

## Cleaner Accessories

<b>ZCD100KIT</b>	Filter and Check Valve Kit
<b>CAP1016S</b>	Quick Release Connector
<b>ZCD100A-DS</b>	Heated Tube Adapter for DreamStation/PR System One
<b>ZCD100A-A10</b>	Heated Tube Adapter Lid for ResMed AirSense 10
<b>ZCD100A-S9</b>	Heated Tube Adapter Lid for ResMed S9
<b>ZCD100TC</b>	Replacement Tube Connector (used when cleaning machine)
<b>ZCD100A-P</b>	Replacement Ozone Plug Adapter (used when not cleaning machine)
<b>ZCD100IT</b>	Replacement Injection Tube with Check Valve
<b>ZCD100F</b>	Replacement Filter
<b>ZCD100PA</b>	Replacement Power Adapter

# CPAP BATTERIES

- Outlet functions like a wall socket
- Universal 110-240V recharging conversion
- 1-Year Warranty

Please note,  
these CPAP  
Batteries have  
a 110V outlet.

## Zopec Explore 5700

Item#  
**CAP1024-2157**

### Specs

- 1 USB-C (5V/3A) and 2 USB (5V/2.1A) ports
- 2.6 lbs
- 8.5 x 6 x 1 in
- 54,200 mAh /160 Wh
- FAA/TSA compliant (carry-on only)



### Zopec Explore Batteries Include

- Universal Wall Charger
- Fast Car Charger
- Travel Pouch
- Run any electronics with power requirement of 120 watts or under
- AutoSwitch and Bypass Technology for uninterrupted sleep when power outages occur
- Unlimited linking of multiple batteries together for longer use before recharging

## Zopec Explore 8200

Item#  
**CAP1024-2182**

### Specs

- 1 USB-C (5V/3A) and 2 USB (5V/2.1A) ports
- 3.5 lbs
- 7.75 x 5.5 x 1.75 in
- 80,000 mAh /288 Wh
- Not FAA/TSA compliant (exceeds power capacity)



## Zopec Voyage

Item#  
**CAP1024-21ZDP02**

Includes battery, wall charger, and travel bag

### Specs

- 1 USB-C (45W) and 2 USB (5V/2.1A) ports
- 2.15 lbs
- 7.25 x 5 x 1.75 in
- 40,000 mAh /146 Wh
- FAA/TSA compliant



## Portable Outlet UPS CPAP Battery

Item#  
**CAP10082**

Includes battery, wall charger, and travel bag

- Autoswitch technology: When power goes out, your CPAP stays on
- Three-prong AC outlet
- Two USB ports
- Digital display
- Convenient travel bag
- Meets TSA/FAA battery guidelines
- 1-year warranty
- Rated for over 500 charge cycles
- Recharges in less than 4 hours
- May be used for other devices powered by AC or USB
- AC and USB ports may be used simultaneously



### Specs

- 2.875 lbs
- 5.375 in wide x 9 in long x 1.5 in tall
- Rated Capacity: 159Wh/14.4Ah
- Rated Voltage: 14.5V

- Charge Input: 16.8V/2A
- USB Output: 5V/4.2A
- AC Output: 110v AC 60hz/100w
- AC Adapter Input: 100V-240V/50-60Hz

# Zopec powers more than CPAP

*Reliable and multiuse backup solutions*

ZOPEC EXPLORE 8200 is ideal for:

## Ventilators

- Adds 9 hours of additional use to internal battery time

## Power Wheelchairs

- A perfectly portable, lightweight solution at just 3 lbs

## Stair Lifts

- Runs for 4 hours continuously (standard backup battery provides just 0-12 runs)



ZOPEC EXPLORE 5700 is ideal for:

## Lift Chairs

- Provides 2-3 hours continuous operation. (Standard backup is two 9V batteries and only lifts one time)



**Bypass** keeps Zopec charged while devices are powered normally.

**Autoswitch** instantly changes the device's power source to Zopec in the event of a power outage.



## Snugz CPAP Mask Liners

Item#  
**CAP7001**

Machine-Washable Liners  
for Most Full Face Masks  
**2/pack**

**CAP7002**

Machine-Washable Liners  
for Most Nasal Masks  
**2/pack**



- 2 per pack, 90-day supply
- Machine washable
- Fits MOST nasal or full face CPAP masks (Not recommended for ResMed Mirage Quattro or Respirationics Amara View and Dreamwear)
- Blend of spandex, moisture-wicking polyester and lycra fibers
- Latex-free

Snugz is a machine washable, one-size-fits-most liner for CPAP masks that is designed to help reduce noisy leaks, prevent redness and irritation, and feel great against the skin. The liner helps to grip the skin and keep the mask in place around the nose and mouth, reducing the need for a chinstrap. The improved seal is especially beneficial for high-pressure/high-flow PAP patients and can reduce nightly sleep disturbances and AHI.


## Soft Headgear Strap Covers

Item#  
**CAP2004**

**NEW**



- Reduces strap marks
- No bulky and uncomfortable velcro
- Reduces strap irritation on sensitive skin
- Thick padding for extra comfort
- No need for clips – just slip wraps around headgear straps
- Easy and helpful for patients with arthritis in hands or fingers
- Can be machine washed and dried
- Wraps securely grip headgear and slide to adjust
- One size fits all
- Made from 230g Coral Velvet



## RemZzzs CPAP Mask Liners

Item#

**CAP3001(S,M,L)**

Liners for Full Face ResMed & Respirationics

OEM#	K1-FS	SMALL
	K2-FM	MEDIUM
	K3-FL	LARGE

**CAP3002(S,M,L)**

Liners for Full Face Fisher & Paykel, Hans Rudolph, and Sunset HCS

OEM#	K4-FS	SMALL
	K5-FM	MEDIUM
	K6-FL	LARGE

**CAP3003(S,M,L)**

Liners for Nasal Respirationics

OEM#	K7-NS	SMALL
	K8-NM	MEDIUM
	K9-NL	LARGE

**CAP3004(S,M,L)**

Liners for Nasal ResMed and Fisher & Paykel

OEM#	K10-NS	SMALL
	K11-NM	MEDIUM
	K12-NL	LARGE

**CAP3005L**

Liners for Respirationics FitLife Mask

OEM#	6A-FLK	SMALL AND LARGE
	6B-FLK	EXTRA LARGE

**CAP3006**

Liners for Nasal Pillows

OEM# 7A-NPK

**CAP3007**

Liners for Respirationics DreamWear and Amara View

OEM# 6C-AVK

### 30-Day Supply

Boxes sold individually.

Full face and nasal mask RemZzzs include 30 facial liners and 2 forehead liners.

- Increases compliance by adding a comfortable and absorbent pad to CPAP mask
- Reduces air leaks and protects skin from red marks



# CASH ITEMS

## CPAP Moisture Therapy

Item# **CAP6100L**  
1oz Tube • 18/cas

- Formulated to prevent chafing, irritation, dryness or cracking from using a CPAP mask or other delivery device
- Apply to the face at the beginning of treatment and as needed
- Petroleum-free with aloe vera, vitamins A & E, and argan oil



## Pur-Sleep CPAP Aromatherapy

Item# **CAP5001**  
VaporClear Aromatherapy Starter Kit  
OEM# VAPOR

### Kit Includes:

- 1 aluminum diffusion stand
- 4 reusable diffusion pads
- 1.5ml bottle of Clear (peppermint, lime)
- 1.5ml bottle of Calm (eucalyptus, lavender, rosemary)



### 30ml Aromatic Refills

Scent/OEM#	Item#
Black Cherry	<b>CAP5004-BCH</b>
Calm	<b>CAP5004-CLM</b>
Clear	<b>CAP5004-CLR</b>
Dark Chocolate	<b>CAP5004-DCH</b>
Deep	<b>CAP5004-DEE</b>
French Vanilla	<b>CAP5004-FRV</b>
Fresh	<b>CAP5004-FSH</b>
Peace	<b>CAP5004-PCE</b>
Spice	<b>CAP5004-SPC</b>
Strawberry	<b>CAP5004-STR</b>

Item# **CAP5002**  
Reusable Diffusion Pads  
10/pack • OEM# DP10

Item# **CAP5003**  
Replacement Aluminum Diffusion Tray  
Each • OEM# TRAY

- Aromatherapy products specifically for use with CPAP
- Set up the diffuser right below the CPAP filter intake, add a drop of essential oil to the pad and the fragrance is drawn into the CPAP machine



## H2Doze Distilled Water for CPAP


Item# **CAPH20-6**  
16.9oz Bottles  
6/pack

- TSA cleared bottle when accompanied with CPAP machine
- Convenient 16.9 oz bottle perfect for travel
- Enough water for multiple nights of use
- Distilled water designed to protect CPAP equipment
- Allows patient to breathe healthy, safe air throughout the night



## Quick Tube Release Connector

Item# **CAP1016**  
Each



## CPAP Final Bacteria Filter

Item# **BF2200**  
22mm Bacteria Disk Filter  
5/pack



Attach between the machine and tube for additional protection against bacteria and allergens

## Circuit Connectors

Item# **RES013**  
22mm OD x 15mm ID/  
22mm OD  
10/pack



Item# **RES014**  
22mm OD x 22mm OD  
With Lip  
10/pack



## Humidifier Tubing Adapter

Item# **RES016**  
22mm ID x 22mm OD  
with Lip  
10/pack



## Oxygen Bleed-In Adapter

Item# **RES019** 10/pack  
22mm with 6mm  
Rotating Stem

**RES019-CAP** 10/pack  
With Stem Cap



## Swivel Connector Silencer

Item# **VEN002**  
For CPAP and Ventilator Circuits



## CPAP Pillows

Item# **CAP4003**

CPAPMax Pillow

Each OEM# 15-501R-1-920



**CPAPMax Replacement Cover**

**CAP4003RC**

Each OEM# 1-630-101-1-933

**Replacement Pillow Case**

**CAP4003C**

Each OEM# 1-626-900-1-02

Item# **CAP4001**

CPAP Pillow

Each

OEM# 14-101R-1-730



**CPAP Pillow Replacement Cover**

**CAP4001RC**

Each OEM# 1-626-100-1-00

**Replacement Pillow Case**

**CAP4001C**

Each OEM# 1-626-800R-1-00




**Cool 3D Mesh Air-Flow**

Reduces heat and perspiration



**Fiberfill**

Soft microfiber pillow top with hypoallergenic fiberfill



**Charcoal Infused High Resiliency Foam**

Charcoal neutralizes odors for freshness every night

### BEFORE

Traditional pillows interfere with mask performance and can cause a stiff neck, shifting, leaks, and facial discomfort.



### AFTER

The CPAPMax Pillow offers superb mask accommodation, excellent support, and unparalleled comfort.



## CPAP Bedside Storage Bag

Item# **CAP1015**



- Well-made and stylish
- Designed to hold CPAP supplies such as a mask and tube
- Weighted lid allows bag to hang off top of nightstand or tabletop



## CPAP Holder

Item# **CAP1010**



- Sits between mattress and boxspring
- Saves space by allowing machine to rest in the area between bed and nightstand
- CPAP machine is kept close to patient's head, allowing more freedom to move throughout the night
- Keeps machine below mask, to help prevent rainout

# CASH ITEMS

## Handheld CPAP Pressure Meter

Item#

**CAP003**

- Measures flow 0-18 cm/H<sub>2</sub>O
- Connects directly to the CPAP tubing without an adapter



## Comfort Tubing Wrap

Item#

**CAP2002**

6ft Nylon  
Tubing Wrap

**Each**

Compare to  
ResMed  
Tubing Wrap



- 6ft nylon wrap with fleece lining maximizes heat retention and reduces noise
- Zipper closure for ease of use

## Deluxe Tubing Wrap

Item#

**CAP2001**

6ft Velour  
Tubing Wrap

**Each**

Compare to  
Respironics  
Tubing Wrap



- 6ft plush velour material adds comfort while preventing rainout
- Zipper closure for ease of use

**Both Tubing Wraps come in attractive retail packaging**



## CPAP RECOMMENDED REPLACEMENT SCHEDULE

Item	Schedule
Heated Circuit Tubing	1 per 3 months
Combination Oral/Nasal Mask	1 per 3 months
Replacement Oral Cushion for Combination Oral/Nasal Mask	2 per month
Replacement Nasal Pillows for Combination Oral/Nasal Mask — Pair	2 per month
Full Face Mask	1 per 3 months
Replacement Full Face Mask Cushion or Liner	1 per month
Replacement Nasal Cushion or Liner	2 per month
Replacement Pillows — Pair	2 per month
Nasal or Pillows Mask	1 per 3 months
Headgear	1 per 6 months
Chinstrap	1 per 6 months
Tubing	1 per 3 months
Filter — Disposable	2 per month
Filter — Nondisposable	1 per 6 months
Replacement Water Chamber	1 per 6 months



## Oxygen

	Page
Pulse Oximeter	21
Firesafe	22
Masks	23
Cannula	24-25
Tubing	25
Cylinder Carriers	26-27
Connectors and Adapters	28
Accessories	29-30
Cleaners	30
Bacteria and HEPA Filters	31
Concentrator Cabinet Filters	32-33
Caire/AirSep/SeQual	32
DeVilbiss/Drive	32
Respironics	33
Nidek	33
Jumao	33

# OXYGEN PRODUCTS

## Finger Pulse Oximeter

Item# **RES5102**

- Oversized 1.3" OLED display face is the largest of any finger pulse oximeter
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO<sub>2</sub>, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- CE, ISO and FDA approved
- Includes lanyard and manual

**NEW**



## Firesafe Thermal Shut-Off Valve

Item#  
**RES010**  
Each


Bidirectional



**NEW**  
**IMPROVED**

- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years

- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

### Fatalities per 100,000 patients



\* 1 death was reported among 73% of the patient population (2013-17)

## FireSafe Kits

Item#  
**RES010KIT**  
2 Firesafe Valves  
1 Tubing Connector

Item#  
**RES010CANNULAKIT**  
2 Firesafe Valves  
1 Tubing Connector  
1 Adult Soft Cannula



## FireSafe Tubing Connector

Item#  
**RES010TC** 100/case  
**RES010TC-10** 10/case



## FireSafe Nozzle

Item#  
**RES010FND-20**  
20/pack



# OXYGEN MASKS

## Simple O<sub>2</sub> Mask

With 7ft Safety Tubing

Item# **RES2100**

Adult Size

50/case

Item# **RES2200V**

Pediatric Size

50/case



HCPCS# A4620

- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## Partial Rebreather O<sub>2</sub> Mask

With 7ft  
Safety Tubing  
& Reservoir Bag

- 5-15 LPM
- High concentration

Item# **RES2101V**

Adult Size

50/case

HCPCS# A4620



## Non-Rebreather O<sub>2</sub> Mask

With 7ft Safety Tubing  
& Reservoir Bag

- 10+ LPM
- Low-resistance  
check valves prevent  
rebreathing and allow  
exhaled gas to escape

Item# **RES2102**

Adult Size

50/case

Item# **RES2202V**

Pediatric Size

50/case



HCPCS# A4620



# OXYGEN CANNULA

## Cannula with Tubing

### Adult Cannula with Tubing

Length	Material	Item#	Qty
4ft	Standard	<b>RES1104</b>	<b>50/case</b>
4ft	Soft	<b>RES1104S</b>	<b>50/case</b>
7ft	Standard	<b>RES1107</b>	<b>50/case</b>
7ft	Soft	<b>RES1107S</b>	<b>50/case</b>
15ft	Standard	<b>RES1115</b>	<b>50/case</b>
25ft	Standard	<b>RES1125</b>	<b>25/case</b>
25ft	Soft	<b>RES1125S</b>	<b>25/case</b>

HCPCS# A4615



- Clean, fresh look
- 1-6 LPM
- Curved nasal prongs for added comfort
- 6-channel safety tubing designed to prevent kinks
- Soft, lightweight material is memory free for long-term comfort

Item#

**RES1207V**

Pediatric Cannula  
with 7ft Supply Tube  
**50/case**



**RES1307V**

Infant Cannula  
with 7ft Supply Tube  
**50/case**



## High Flow Soft Cannula

Length	Item#	Qty
4ft	<b>RES1104SHF</b>	<b>25/case</b>
7ft	<b>RES1107SHF</b>	<b>25/case</b>
15ft	<b>RES1115SHF</b>	<b>25/case</b>
25ft	<b>RES1125SHF</b>	<b>25/case</b>

HCPCS# A4615




- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

## Cannula without Tubing


### Adult Cannula without Tubing

Port	Material	Item#	Qty
Male	Standard	<b>RES1100</b>	<b>50/case</b>
Male	Soft	<b>RES1100S</b>	<b>50/case</b>
Female	Soft	<b>RES1100FS</b>	<b>50/case</b>

HCPCS# A4615



Male port



Female port

## Soft Cannula with Universal Connector

Item# **RES1107SUC**

Adult Soft Cannula with  
Universal Connector  
and 7ft Supply Tube  
**50/case**

HCPCS# A4615




**NEW**

## Dual-Lumen (Demand) Cannula

Length	Item#	Qty
4ft	<b>RES4144V</b>	<b>25/case</b>
7ft	<b>RES4147V</b>	<b>25/case</b>

HCPCS# A4615



- 1-6 LPM
- Designed for use with dual port oxygen conservation delivery systems
- Dual tubing set provides independent tubes for oxygen delivery and sensing of patient breathing
- Smoothbore tubing

## Tabbed Cannula

Item# **RES1007V**

Adult Tabbed Cannula  
with 7ft Supply Tube  
**50/case**

HCPCS# A4615






# OXYGEN CANNULA AND TUBING

## Humidoflo High Flow Nasal Cannula

10/case

Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>

**NEW**




- Includes standard ISO connector and F&P adaptor so it can be used with almost any machine on the market
- Kink-free prongs and flexible, low-noise tubing
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)
- Can be used with F&P 950, 850, or Airvo system

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O2 concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

## Secure Cannula Stickers

Item# **RES080**  
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



## Comfort Cannula with Ear Cushions

Item# **RES1107EC**  
7ft Cannula with  
Ear Cushions  
**25/case**  
HCPCS# A4615



- 1-6 LPM
- Soft foam cushions add comfort and increase patient satisfaction
- Cushions are added during the manufacturing process so they stay in place and no assembly is required

Item# **RES020**  
Cushions Only  
**50 Pairs/case**

- Add cushions to existing cannula



# OXYGEN TUBING

## Oxygen Supply Tube



Length	Item#	Qty
7ft	<b>RES3007</b>	<b>50/case</b>
15ft	<b>RES3015</b>	<b>50/case</b>
25ft	<b>RES3025</b>	<b>25/case</b>
40ft	<b>RES3040</b>	<b>20/case</b>
50ft	<b>RES3050</b>	<b>20/case</b>

HCPCS# A4616

- Crush resistant with low memory to prevent kinking
- Great durability and comfort level
- Latex-free
- 1-6 LPM
- 6 channels

## Oxygen Supply Tube — Green



Length	Item#	Qty
7ft	<b>RES3007G</b>	<b>50/case</b>
15ft	<b>RES3015G</b>	<b>50/case</b>
25ft	<b>RES3025G</b>	<b>25/case</b>
40ft	<b>RES3040G</b>	<b>20/case</b>
50ft	<b>RES3050G</b>	<b>20/case</b>

HCPCS# A4616



# CYLINDER BAGS

## Comfort Shoulder Bags

Cylinder	Item#
ML6/M7	<b>CSBML6M7</b>
C/M9	<b>CSBM9</b>
B/M6	<b>CSBM6</b>
D	<b>CSBD</b>



**CASH ITEM**

- Comfort flex shoulder strap
- Double zipper for easy access
- Machine washable



CSBML6M7



CSBM9



CSBM6



CSBD

## Flex 02 Bag

Item# **CSBM6Flex**  
Reconfigurable Bag  
for M6 and M9  
Cylinders  
Each

- Adjusts to 5 carrying modes



**CASH ITEM**

## Backpack

Cylinder	Item#
B/M6/C/M9	<b>CSBM6BP</b>
D	<b>CSBDBP</b>



**CASH ITEM**



- Fully adjustable and comfortable backpack straps
- Machine washable
- Converts to either backpack or shoulder bag

## Wheelchair/Scooter Carrier

Item# **CSBE**  
For E Cylinders  
Each



**CASH ITEM**



- Multi-part stabilizer straps for easy chair attachment
- Machine washable

# CYLINDER CARRIERS

## Single Tank Rolling Bag/Backpack

Item#  
**CSBD1Roller**  
For D/E Cylinders  
Each



**CASH ITEM**



## Double Tank Rolling/Shoulder Bag

Item#  
**CSBD2Roller**  
For D/E Cylinders  
Each



**CASH ITEM**

- May also be used with shoulder strap



## ML6 Hardshell Tote

Item#  
**CSBML6-Hardshell**  
For Homefill Cylinders  
Each




**CASH ITEM**



## Cylinder Rack

Item# **RES038B**  
Each  
For D, E, or M9 Cylinders —  
Lightweight, Durable  
Plastic Composite Rack



## Oxygen Cylinder Carrier

Item#  
**CSBM6-Carrier**  
Each  
Soft Polyester Tote  
for M6 Cylinders



- Tote handles
- Padded shoulder strap
- Fastening strap for mobility
- Holds up to 8 M6 tanks

## D/E Cart

Item#  
**RES060**  
4/case

- Accommodates one D or E tank
- Sturdy yet lightweight (4.18 LBS welded metal construction)
- Adjustable height handle with PVC grip
- Large EVA Foam wheels have just the right amount of resilience for good traction and maneuverability
- Elegant, grey satin finish
- Bottom fork rest with rubber caps



# OXYGEN CONNECTORS AND ACCESSORIES

## Christmas Tree Connector

Item# **RES001**  
Green without Swivel  
**10/pack**



Item# **RES002**  
Green with Swivel  
**10/pack**



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

## Oxygen Tubing Y-Connector

Item# **RES017**  
Fits 3/16in and 1/4in ID Tubing  
(All 3 sides are the same)  
**10/pack**



- Connects 2 oxygen sources for higher flows

## Tubing Reducer Connector

Item# **RES011**  
Fits 3/16in, 1/4in,  
5/16in, 3/8in, and 7/16in  
ID Tubing  
**10/pack**



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

## Oxygen In Use Signs

Item# **RES003**  
6in x 9in  
With String  
**50/pack**



- Double-sided (English and Spanish)

## Oxygen Accessory Kit

Item# **RES1000**  
RES002 Christmas Tree Connector - Swivel  
RES012 Tubing Connector - Non-Swivel  
RES017 Tubing Y-Connector  
RES018 Tubing Connector - Swivel  
RES032 Wrench  
RES035 Brass Washer

Each



**CUSTOMIZABLE KIT**



- Order all common oxygen accessories in one easy-to-handle kit
- Contents can be customized to fit your needs

## Oxygen Tubing Connector — Non-Swivel

Item# **RES012**  
Fits 3/16in and 1/4in  
ID Tubing  
**10/pack**




- Barbed ends ensure secure, leak-free connections

## Oxygen Tubing Connector — Swivel

Item# **RES018**  
Male/Male **10/pack**

**RES018-25** 25/pack

**RES018KIT-1**  
Female/Male **Each**



- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing

## Washers with Rubber Ring

Item# **RES035**  
Brass  
**25/pack**



Item# **RES036**  
Aluminum  
**25/pack**





# OXYGEN ACCESSORIES

## O<sub>2</sub> Starter Kit

### Item# **RES02KIT**

RES1107 Cannula with 7ft Tube x2  
 RES3015 15ft Supply Tube x2  
 RES002 Christmas Tree Connector — Swivel  
 RES003 Oxygen In Use Sign  
 RES012 Tubing Connector — Non-Swivel  
 RES017 Tubing Y-Connector  
 RES018 Tubing Connector — Swivel  
 RES032 Wrench  
 RES035 Brass Washer



**CUSTOMIZABLE KIT**

## Disposable Bubble Humidifier 500 ml

### Item#

### **RES050-50**

Up to 15 LPM • 6 psi • 500 ml

**50/case**

HCPSC# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

## Wrenches

### Item# **RES031**

With Chain **10/pack**

### Item# **RES032**

Without Chain **10/pack**

### Item# **RES031M**

Metal Wrench with Chain for Post Valve Oxygen Cylinders and CGA 540 Style Regulators **5/pack**



- For use with post valve oxygen cylinders: A, B, C, D, JD or E tanks
- Made of tough, stress-resistant nylon

## Humidifier Connectors

### Item# **HUM008**


L-Connector 10in • **10/pk**

### Item# **HUM008L**

L-Connector 21in • **10/pk**

### Item# **HUM008S**

Straight Connector 9in • **10/pk**



## Water Traps

### Item# **RES037**

Push On End Caps • **Each**

### Item# **RES037T**

Threaded End Caps • **Each**



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

# ACCESSORIES AND CLEANERS

## Homefill Tubing Connector


Item# **RES039**  
10/pack

- Compatible with Invacare, Drive, and Respiration homefill machines
- Durable, coiled tubing




## Homefill Racks

Item# **RES039-RACKB**  
For Invacare Perfecto2  
4/case



**RES039-RACKW**  
For Invacare Platinum  
4/case

**RES039-RACKD**  
For DeVilbiss 525/1025 Series  
4/case



Easily attaches onto concentrator and is designed to conveniently hold Homefill and transfill oxygen systems.

## Oxygen Liter Meter

Item# **RES034**  
0-10 Liter/Minute  
Each



- Oxygen pen
- Works with tanks and concentrators

## RoEzIt Dermal Care

Item# **RES6104L**  
1oz Tube • 18/cas



**CASH ITEM**

- Apply before beginning oxygen therapy and as needed during treatment
- Moisturizes the nose, lip, and over the ear where friction may cause discomfort from cannula
- Petroleum-free with aloe vera and vitamins A and E

## Germicidal Deodorizing Cleaner

Item# **RES6100**  
12/case  
22oz spray bottles



- Heavy duty cleaner contains no bleach, alcohol, or phenols
- An effective bactericide, virucide, and fungicide
- Formulated to disinfect hard, non-porous surfaces

## Germicidal Deodorizing Cleaner — Gallon

Item# **RES6102**  
1 Gallon




- No-scrub, hospital strength germicidal cleaner
- Cleans, deodorizes, and disinfects surfaces, equipment and non-critical instruments

# OXYGEN BACTERIA AND HEPA FILTERS

## EverFlo Intake Filter

Item# **BF1001**  
Each


• Respiration EverFlo



## Extended Life Bacteria Filter with Muffler

Item# **OF6011**  
Each

• DeVilbiss/Drive 525 DS Series



## Extended Life Bacteria Filter

Item# **BF5000**  
Each

• Nidek Nuvo



## Extended Life Intake Bacteria Filter

Item# **OF1025D**  
Each

• DeVilbiss/Drive 1025



## Extended Life Intake Bacteria Filter

Item# **BF1103**  
Each

• CAIRE Companion 5



## Intake Bacteria Filter

Item# **BF7000**  
Each

• Jumao



## Final Bacteria Filter


Item# **BF516**  
**10/pack** 1/4in Barbed on Both Sides

• Sequal Integra



## Final Bacteria Filter — 1/4in to 3/8in

1/4in —  
3/8in —



Item# **BF1438**  
**10/pack**  
1/4in to 3/8in Stepped Barb  
on Both Sides

• AirSep, DeVilbiss/Drive, and Nidek units

## Final Bacteria Filter — 3/16in




Item# **BF316**  
**10/pack**  
3/16 Stepped Barb and 1/4in  
Straight Barb on Both Sides

• Respiration and 3B units

# CAIRE/AIRSEP/SEQUAL


## CAIRE Companion 5



Item# **BF1103**  
Extended Life Intake  
Bacteria Filter  
**Each**  
OEM# 2185


**NEW**

## SeQual Eclipse



Item# **OF11006**  
Foam Cabinet Filter  
**10/pack**  
21/2" x 21/2" x 3/4"

## AirSep Newlife, QuietLife, Intensity



Item# **OF8002**  
Foam Cabinet Filter  
**10/pack**  
51/4" x 47/8" x 7/8"  
OEM# F1-50-08-05



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in  
Barbed  
**10/pack**  
OEM# F1-30-03



Item# **OF8001**  
Felt Filter  
**10/pack**  
17/8" OD x 3/4" ID x 1"  
OEM# M1 161-1



**KIT**

Item# **OF8002SP**  
Newlife Kit  
2 - OF8002  
1 - BF1438  
1 - OF8001  
Each

# DeVILBISS/DRIVE

## 1025 Series



Item# **OF1025D**  
Extended Life Bacteria  
Filter  
**Each**



Item# **OF5004**  
Foam Cabinet Filter  
**10/pack**  
33/4" x 33/4" x 1/2"  
OEM# 303DZ-605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651


## 525 DS Series



Item# **OF5004**  
Foam Cabinet Filter  
**10/pack**  
33/4" x 33/4" x 1/2"  
OEM# 303DZ-605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



Item# **OF6011**  
Extended Life Bacteria  
Filter — 22mm fittings  
**Each**



# RESPIRONICS

## EverFlo



Item# **OF1003SP**  
EverFlo Kit  
1 - BF1001  
1 - BF316  
Each



KIT



Item# **OF1006SP**  
Cabinet Filter Kit  
(No Cover)  
**2 Filters, 1 Casing**  
OEM# 1055625



Item# **OF1006**  
Cabinet Filter  
**10/pk**  
5 11/16" X 4 1/2" X 7/16"



Item# **OF1006-CASE**  
Foam Cabinet Filter  
Casing  
**10/pack**



Item# **BF1001**  
EverFlo Intake Filter  
**Each**  
OEM# 1038831



Item# **HUM008**  
Humidifier  
L-Connector  
With 10in  
Tube  
**10/pk**



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and  
1/4in Straight Barb  
**10/pack**  
OEM# H621



Item# **OF2000**  
Sieve Bed Filter 4.125in  
**5/pack**

# NIDEK

## Nuvo



Item# **OF11003SP**  
Nuvo Kit  
2 - OF11003  
1 - BF1438  
1 - BF5000  
Each



KIT



Item# **OF11003**  
Foam Cabinet Filter  
**10/pack**  
9 1/2" X 4" X 1/4"  
OEM# 9250-1025



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# 7631-1053



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**  
OEM# 9250-1028

# JUMAO

## Jumao



Item# **OF7000**  
Foam Cabinet Filter  
**10/pack**



Item# **BF7000**  
Bacteria Filter  
**Each**

## Other

## Page

Aerosol	35
Handheld Nebulizers	36
Nebulizer Kits	37
Nebulizer and Anesthesia Filters	38
Suction	39-40
Closed Suction	41-42
Tracheostomy	43
Ventilator Filters, Masks, and Connectors	44
Ventilator Circuits	45-46
StrongArm Comfort Cane	47
Anesthesia	48-50
Anestand	48-49
GuideBlade	50

# OTHER

# AEROSOL

## Aerosol Mask

Item# **RES2110**  
Adult Size  
**50/case**

Item# **RES2210**  
Pediatric Size  
**50/case**

HCPCS# A7015



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



## Peak Flow Meter

Item# **RES052**  
**10/case**

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**  
Pediatric Mouthpiece  
for Peak Flow Meter  
**10/case**

## Drainage Bag with T-Adapter

Item# **RES070**  
750cc • **50/case**



## Drainage Bag with Y-Adapter

Item# **RES071**  
750cc • **50/case**



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

## SperaChamber

Anti-Static Valved Inhaler Chamber  
for Metered Dose Inhaler

Chamber Only

Item#  
**RES195-10 10/pack**  
**RES195-50 50/pack**

Chamber with Mask  
**10/pack**

Item#  
**RES195MS Small**  
**RES195MM Medium**

**RES195ML Large (Adult)**

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free
- Dog mask appeals to kids



# HANDHELD NEBULIZERS

## Handheld Mesh Nebulizer

Item#  
**NEB402**

**NEW**

**With internal rechargeable battery**

**Included**

Medication cup  
USB cord  
Mask  
Mouthpiece  
Portable bag  
Manual



**Particle Size** MMAD < 5µm  
**Nebulization Rate** 0.25ml/min-0.9ml/min  
**Medication Cup Capacity** 10ml (Max)  
**Device Weight** Approx. 1/4 lb  
**Device Size** Approx. 1 1/2 x 2 x 3 7/8 inches

## Mini Mesh Nebulizer

Item#  
**NEB400**

**Included**

USB cord  
Medicine cup  
Mouthpiece  
Adult mask  
Pediatric mask  
Carrying bag

- Runs on USB or 2 AA batteries
- Quiet
- Fast treatment times
- Small particle size
- Attainable for more patients

**Replacement Parts**

<b>NEB400CUP</b>	Medication Cup
<b>NEB400MP</b>	Mouthpiece
<b>NEB400MASKP</b>	Pediatric Mask
<b>NEB400MASKA</b>	Adult Mask



## Handheld Compressor Nebulizer

Item#  
**NEB200**

**Included**

Compressor  
Medication cup  
AC power adapter  
Air hose  
Elbow  
Mouthpiece  
Mask  
4 AA batteries



- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required




# NEBULIZER KITS

## Disposable Nebulizer Kit with T-Piece

Item# **RES091**  
Disposable  
Single Patient  
**50/case**  
HCPCS# A7003

- Includes:
- 6cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube




- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

## Reusable Nebulizer Kit with T-Piece

Item# **RES093**  
Reusable — Single Patient  
**10/case**  
HCPCS# A7005

- Includes:
- 10cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

## Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**  
Disposable  
Single Patient  
**50/case**  
HCPCS# A7003

- Includes:
- 6cc jet nebulizer bottle
  - Angled mouthpiece
  - 7ft supply tube



## Disposable Nebulizer Kit with Mask

Item# **RES090**  
Adult Size  
**50/case**

Item# **RES092**  
Pediatric Size  
**50/case**

HCPCS# A7003

### Includes

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

Disposable single patient use for long term delivery procedures

6cc bottle is easy to assemble and disassemble



# NEBULIZER FILTERS

## Plastic Intake Filter

Item# **BF015**

**10/pack**

HCPCS# A7013



## Plastic Intake Filter with Head

Item# **BF016**

**10/pack**

HCPCS# A7013



## Plastic Intake Filter with Head

Item# **BF017**

**5/pack**

HCPCS# A7013



## Foam Filter

Item# **BF018**

**10/pack**

HCPCS# A7013



## Plastic Intake Filter for Respironics InnoSpire

Item# **BF026S**

**10/pack**

HCPCS# A7013

- Mini
- MicroElite
- MiniElite



## Nebulizer Filter Compatibility

### Caire

Medi-Mist	BF018
-----------	-------

### DeVilbiss

PulmoNeb	BF017
PulmoNeb LT	BF017
Sunrise Compact	BF018
Pulmo-Aide	BF017
Pulmo-Aide Compact	BF017
PulmoMate	BF017

### Evo

AeroMax	BF018
SportsNeb	BF018
Super SportsNeb	BF018

### Invacare

Passport	BF018
Envoy	BF018
Envoy Jr.	BF018
Select	BF018

### John Bunn

Nebulite EV	BF018
LT2	BF018
LX	BF018

### MedLine

SportsMist	BF018
------------	-------

### Respironics

Inspiration 323	BF015
Inspiration 626 & 929	BF016
MisterNeb	BF016
Inspiration Elite	BF016
MiniElite	BF018
InnoSpire	BF026S

### Salter

AIRE Plus	BF018
-----------	-------

### Schuco

2000	BF018
3000	BF018

### Sunset

NEB100 Series	BF018
NEB300 Series	BF018

# SUCTION

## Suction Machine

Item#

**SU100DC**

With Internal  
Rechargeable  
Battery

**SU100AC**

Without Battery



### Comes with

Machine  
Carry Bag  
800cc Suction Canister  
AC Adapter & Power Cord  
6ft Suction Tubing  
3.5in Suction Tube Connector  
4.5in Suction Tube Connector  
BF1438S Bacteria Filter

- Effective removal of respiratory fluids
- Durable yet lightweight
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty

### Accessories and replacement parts

**SU100DC-PA** DC Mobile Power Adapter (for car power outlet)

**SU100DC-Adapt** Replacement Power Adapter (without cord)

**SU100DC-Cord** Replacement Power Cord (without adapter)

**SU100DC-BAG** Carrying/Storage Bag

**SU100-PA** Car Charger (DC Mobile Power Adapter)

**SU100DC-BATTERY** Internal rechargeable battery replacement

## Suction Kit for SU100DC

Item# **RES026S-SS**

RES023AW Suction Canister  
RES024 Suction Tubing Connector 3.5in  
RES024M Suction Tubing Connector 4.5in  
BF1438S Bacteria Filter  
RES025 Suction Tubing 1/4in x 6ft

HCPCS#

CANISTER - A7000

TUBING - A7002



COMPATIBLE WITH  
DRIVE  
DEVILBISS  
& SUNSET  
MACHINES

## Suction Filter Component Kit for SU100DC

Item# **RES024-MKIT**

RES024 Suction Tubing Connector 3.5in  
RES024M Suction Tubing Connector 4.5in  
BF1438S Bacteria Filter

COMPATIBLE WITH  
DRIVE  
DEVILBISS  
& SUNSET  
MACHINES

## Suction Filter Component Kit

Item# **RES024XL-3KIT**

RES024XL Suction Tubing Connector 16in  
BF3812 Filter  
RES028 Elbow Connector

COMPATIBLE WITH  
MEDLINE  
& PRECISION  
MACHINES

## Suction Canister

Item#

**RES023AW**

800cc • Single Use  
**48/case**

HCPCS# A7000

Item#

**RES025AW1200**

1200cc • Single Use  
**48/case**

COMPATIBLE WITH  
DRIVE  
DEVILBISS  
& SUNSET  
MACHINES



- Float valve shutoff

## Suction Kit

Item# **RES026S-E**

RES023AW Suction Canister  
RES024M Suction Tubing Connector 4.5in  
RES025 Suction Tubing 6ft  
RES028 Elbow Connector  
BF3812 Bacteria Filter

HCPCS#

CANISTER - A7000

TUBING - A7002



COMPATIBLE WITH  
DRIVE  
DEVILBISS  
& SUNSET  
MACHINES

# SUCTION

## Suction Machine Inline Filters

Item# **BF1438S**  
1/4in to 3/8in  
10/pack



Item# **BF3812**  
3/8in to 1/2in  
10/pack



Item# **BF1312**  
1/3in to 1/2in  
1/pack



Item# **BF516TS**  
5/16in Barbed with  
Threaded NPT  
10/pack



## Suction Tubing

### 1/4in

Length	Item#	Qty
18in	<b>RES025S</b>	<b>50/case</b>
6ft	<b>RES025</b>	<b>10/case</b>

### 3/16in

Length	Item#	Qty
18in	<b>RES025-31618</b>	<b>100/case</b>
20in	<b>RES025-31620</b>	<b>100/case</b>
6ft	<b>RES025-31672</b>	<b>50/case</b>

HCPCS# A7002



## Suction Tube Connectors

### 10/pack

Size	Item#
1.5in	<b>RES024S</b>
3.5in	<b>RES024</b>
4.5in	<b>RES024M</b>
6in	<b>RES024L</b>
10in	<b>RES024L10</b>
16in	<b>RES024XL</b>

**NEW**



## Rubber Suction Elbow Connector Drive DeVilbiss Style

Item#  
**RES028**  
6/pack



## Suction Elbow Connector Precision Medical Style

Item#  
**RES033**  
6/pack



## Bulb Tip Yankauers

Item# **TRAA001S**  
Vented • **50/case**

Item# **TRAA002S**  
Non-Vented • **50/case**

HCPCS# A4628



- Sterile
- Single patient use
- Not made with natural latex rubber



## Standard

### 24 Hour Adult

#### Endotracheal 540mm

Item#	Description
<b>VEN46001000A</b>	10Fr
<b>VEN46001000AM</b>	10Fr, MDI Port
<b>VEN46001200</b>	12Fr
<b>VEN46001200M</b>	12Fr, MDI Port
<b>VEN46001400</b>	14Fr
<b>VEN46001400M</b>	14Fr, MDI Port
<b>VEN46001600</b>	16Fr
<b>VEN46001600M</b>	16Fr, MDI Port

#### Tracheostomy 305mm

<b>VEN46011000A</b>	10Fr
<b>VEN46011000AM</b>	10Fr, MDI Port
<b>VEN46011200</b>	12Fr
<b>VEN46011200M</b>	12Fr, MDI Port
<b>VEN46011400</b>	14Fr
<b>VEN46011400M</b>	14Fr, MDI Port
<b>VEN46011600</b>	16Fr
<b>VEN46011600M</b>	16Fr, MDI Port

### 72 Hour Adult

#### Endotracheal 540mm

Item#	Description
<b>VEN56051000A</b>	10Fr
<b>VEN56051000AM</b>	10Fr, MDI Port
<b>VEN56051200</b>	12Fr
<b>VEN56051200M</b>	12Fr, MDI Port
<b>VEN56051400</b>	14Fr
<b>VEN56051400M</b>	14Fr, MDI Port
<b>VEN56051600</b>	16Fr
<b>VEN56051600M</b>	16Fr, MDI Port

#### Tracheostomy 305mm

<b>VEN56061000A</b>	10Fr
<b>VEN56061000AM</b>	10Fr, MDI Port
<b>VEN56061200</b>	12Fr
<b>VEN56061200M</b>	12Fr, MDI Port
<b>VEN56061400</b>	14Fr
<b>VEN56061400M</b>	14Fr, MDI Port
<b>VEN56061600</b>	16Fr
<b>VEN56061600M</b>	16Fr, MDI Port




## Neonatal/Pediatric With Y-Adapter

Item#	Description
<b>VEN66000500</b>	305mm, 5Fr
<b>VEN66000600</b>	305mm, 6Fr
<b>VEN66000700</b>	305mm, 7Fr
<b>VEN66000800</b>	305mm, 8Fr
<b>VEN66001000</b>	405mm, 10Fr
<b>VEN66000500M</b>	305mm, 5Fr, MDI Port
<b>VEN66000600M</b>	305mm, 6Fr, MDI Port
<b>VEN66000700M</b>	305mm, 7Fr, MDI Port
<b>VEN66000800M</b>	305mm, 8Fr, MDI Port
<b>VEN66001000M</b>	405mm, 10Fr, MDI Port

## Deluxe

### Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single handed operation design, which allows you to keep your eye on the patient and do no harm. One handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.



### Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.

**MDI Port** provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

**Irrigation Port** offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

**Isolation Valve** seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

**Black Ring Indicator** ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.

# CLOSED SUCTION

## Deluxe

### 24 Hour Adult

#### Endotracheal 540mm

Item#	Description
<b>VEN160100A</b>	10Fr
<b>VEN160100AM</b>	10Fr, MDI Port
<b>VEN160120</b>	12Fr
<b>VEN160120M</b>	12Fr, MDI Port
<b>VEN160140</b>	14Fr
<b>VEN160140M</b>	14Fr, MDI Port
<b>VEN160160</b>	16Fr
<b>VEN160160M</b>	16Fr, MDI Port

#### With Wedge + Adapter

Item#	Description
<b>VEN162100A</b>	10Fr
<b>VEN162100AM</b>	10Fr, MDI Port
<b>VEN162120</b>	12Fr
<b>VEN162120M</b>	12Fr, MDI Port
<b>VEN162140</b>	14Fr
<b>VEN162140M</b>	14Fr, MDI Port
<b>VEN162160</b>	16Fr
<b>VEN162160M</b>	16Fr, MDI Port

#### Tracheostomy 305mm

<b>VEN161100A</b>	10Fr
<b>VEN161100AM</b>	10Fr, MDI Port
<b>VEN161120</b>	12Fr
<b>VEN161140</b>	14Fr
<b>VEN161160</b>	16Fr
<b>VEN161160M</b>	16Fr, MDI Port

#### With Wedge + Adapter

<b>VEN163100A</b>	10Fr
<b>VEN163100AM</b>	10Fr, MDI Port
<b>VEN163120</b>	12Fr
<b>VEN163140</b>	14Fr
<b>VEN163160</b>	16Fr

### 72 Hour Adult

#### Endotracheal 540mm

Item#	Description
<b>VEN365100A</b>	10Fr
<b>VEN365100AM</b>	10Fr, MDI Port
<b>VEN365120</b>	12Fr
<b>VEN365120M</b>	12Fr, MDI Port
<b>VEN365140</b>	14Fr
<b>VEN365140M</b>	14Fr, MDI Port
<b>VEN365160</b>	16Fr
<b>VEN365160M</b>	16Fr, MDI Port

#### With Wedge + Adapter

Item#	Description
<b>VEN367100A</b>	10Fr
<b>VEN367100AM</b>	10Fr, MDI Port
<b>VEN367120</b>	12Fr
<b>VEN367120M</b>	12Fr, MDI Port
<b>VEN367140</b>	14Fr
<b>VEN367140M</b>	14Fr, MDI Port
<b>VEN367160</b>	16Fr
<b>VEN367160M</b>	16Fr, MDI Port

#### Tracheostomy 305mm

<b>VEN366100A</b>	10Fr
<b>VEN366100AM</b>	10Fr, MDI Port
<b>VEN366120</b>	12Fr
<b>VEN366120M</b>	12Fr, MDI Port
<b>VEN366140</b>	14Fr
<b>VEN366140M</b>	14Fr, MDI Port
<b>VEN366160</b>	16Fr
<b>VEN366160M</b>	16Fr, MDI Port

#### With Wedge + Adapter

<b>VEN368100A</b>	10Fr
<b>VEN368100AM</b>	10Fr, MDI Port
<b>VEN368120</b>	12Fr
<b>VEN368120M</b>	12Fr, MDI Port
<b>VEN368140</b>	14Fr
<b>VEN368140M</b>	14Fr, MDI Port
<b>VEN368160</b>	16Fr
<b>VEN368160M</b>	16Fr, MDI Port

### Neonatal/Pediatric

Item#	Description
<b>VEN26000500</b>	305mm, 5Fr
<b>VEN26000500M</b>	305mm, 5Fr, MDI Port
<b>VEN26000600</b>	305mm, 6Fr
<b>VEN26000600M</b>	305mm, 6Fr, MDI Port
<b>VEN26000700</b>	305mm, 7Fr
<b>VEN26000700M</b>	305mm, 7Fr, MDI Port
<b>VEN26000800</b>	305mm, 8Fr
<b>VEN26000800M</b>	305mm, 8Fr, MDI Port
<b>VEN26001000</b>	405mm, 10Fr
<b>VEN26001000M</b>	405mm, 10Fr, MDI Port

# TRACHEOSTOMY

## Tracheostomy Mask

Item#

**RES2130**

Adult Size

**50/case**

**RES2230**

Pediatric Size

**50/case**

HCPCS# A7525



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

## HMEs

### Heat Moisture Exchangers



**RES027FLEX**

Flex HME —  
Foam



**RES027**

T-HME



**RES027OP**

T-HME with  
Oxygen Port



**RES027IL**

Inline HME —  
Large Volume



**RES027NE**

Inline HME —  
Neonatal

TV Range (mL)	300-1500	250-1000	250-1000	15-50	15-50
Moisture Output (mg H <sub>2</sub> O/L)	33.8, Vt = 0.5L	27.3, Vt = 0.5L	26, Vt = 0.5L	31, Vt = 0.1L	30, Vt = 0.2L
Weight (g)	17	5	6.5	10	5
Resistance (cm H <sub>2</sub> O)	0.5@30 LPM, 1.34@60 LPM	0.7@60 LPM	0.4@30 LPM	0.4@10 LPM	0.5@10 LPM
Dead Space (mL)	90	12.5	12	13	2.5
Connection (mm)	22M, 15F	15F	15mmF, 5.5mm Oxygen port	15M, 22M/15F	15M, 15F

### Tracheostomy Ties with Fastener Tabs

**20/case**


One-piece

- Easy on

Item#

**RES029** Adult

**RES029P** Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton



Two-piece

- Easy off

Item#

**RES030** Adult

**RES030P** Pediatric



# VENTILATOR FILTERS, MASKS, AND CONNECTORS

## NIV Mask with Headgear — Non-vented

Each

Item#

**CM038M** Medium

**CM038L** Large

- Disposable, single patient use
- PC frame with PVC air cushion seal
- Low dead space for low CO2 buildup
- Clear mask allows a good view of mouth and nasal area
- Headgear with Velcro fasteners provide maximum versatility and movement
- Dual 360 degree swivels at the mask attachment and at the standard 22mm tubing attachment point



## Bacteria Filter

Item#

**BF2200**

22mm Bacteria  
Disk Filter

**5/pack**



## Respironics Trilogy Ventilator Filter

Item#

**CF1009F**

**2/pack**

OEM# 1029330



## Air Inlet Vent Filter for Philips Respironics V60

Item#

**VF002**

For use with Philips Respironics V60

**5/pack**



## Inlet Vent Filter for Pulmonetic LTV Series

Item#

**VF001**

For use with Pulmonetic LTV Series

**5/pack**



## Valved T-Adapter — Adult

Item#

**VEN001**

22F-18F-22M/15F

**Each**



## Ventilator Felt Filter

Item#

**VF006**

For use with NPB LP6 & LP10

**10/pack**



## Swivel Connector Silencer

Item#

**VEN002**

For Ventilator Circuits and CPAP

**Each**



## Circuit Connectors

Item#

**RES013**

15mm ID x 22mm OD

**10/pack**

**RES014**

22mm OD x 22mm OD

**10/pack**



## Fixed Exhalation Port

Item#

**VEN003**

For Ventilator Circuits and CPAP  
22 mm OD x 15 mm ID

**Each**

**NEW**



# VENTILATOR CIRCUITS WITHOUT HUMIDIFIER CHAMBER

## Passive Heated Breathing Circuit


Item#

**VC1114**

For use with MR850 Humidifier

**10/case**

OEM# RT1114



## Passive Non-Heated Breathing Circuits

For use with Respiration Trilogy

Item#

**VC1000**

With swivel exhalation port

**Each**

OEM# 1073228 & 1073236



**VC1100**

With fixed exhalation port

**Each**

OEM# 1069467 & 1069468

## Cough Assist Kits

Item#

**CAKIT001**

**Each**

OEM# 1090835

**NEW**

- 1 RES013 Circuit Connector  
22mm OD, 15mm ID x  
22mm OD
- 1 BF2200 Bacteria/Viral  
Filter
- 1 TUB06 Circuit Tube
- 

**CAKIT002**

**Each**

- 1 RES013 Circuit Connector  
22mm OD, 15mm ID x  
22mm OD
- 1 BF2200 Bacteria/Viral  
Filter
- 1 TUB06 Circuit Tube
- 1 RES019 Oxygen  
Enrichment Adapter





# VENTILATOR CIRCUITS WITHOUT HUMIDIFIER CHAMBER


## Non-Heated Breathing Circuits

For use with Respiration Trilogy  
22mm cuffs  
**30/case**




Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
<b>VC1001</b>	1073223 & 1073232	Active	No	Adult	22mm
<b>VC1001P</b>	1073220 & 1073229	Active	Yes	Pediatric	15mm
<b>VC1004</b>	1069443 & 1069444	Passive	Yes	Adult	22mm
<b>VC1004P</b>	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#  
**VC1001**  
Active Adult  
OEM#  
1073223  
1073232




Item#  
**VC1001P**  
Active Pediatric  
OEM#  
1073220  
1073229




To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.



Item#  
**VC1004**  
Passive Adult  
OEM#  
1069443  
1069444



Item#  
**VC1004P**  
Passive Pediatric  
OEM#  
1073224  
1073233



# STRONGARM

## StrongArm® Comfort Cane

**NEW**

- Ergonomic handgrip supports your forearm and relieves pressure from your wrist
- Offset bend centers grip over tip for maximum stability when cane is at the user's side
- Lightweight & strong
- 12 height adjustments
- Comfort Fit Handle
- Supports up to 500 lbs
- For left & right arm
- Self-standing base: Larger 'footprint' greatly improves traction and overall stability
- Safely walk on grass, gravel, sand, carpet, hard surfaces, and stairs



**10/case**

**BM001** Black

**BM001-BL** Blue

**BM001-TI** Titanium

**One of each color pack**

**BM001-3MC**



**STRONGARM COMFORT CANE**





## Customized anesthesia workstation for every anesthetizing location

Item#

**ANE77006**

Anestand Anesthesia Stand, Complete Tray and Clamp

**ANE77001** Replacement Tray

**ANE77005** Replacement Clamp

### Benefits to Providers

- Keep medications near IV port
- Ergonomical work space
- Makes wiping IV port easy
- Extra set of hands when help not around
- Helps make RSI safe
- Creates valuable space for non-OR anesthesia

### Value to Hospital

- Helps with regulatory compliance
- Sharps management
- Could help reduce HAIs
- Demonstrates support for anesthesia team
- Reusable/single fixed cost

1. The recessed tray holds and organizes supplies for induction of anesthesia and most anesthetic procedures
2. Clips keep IV tubing securely in place, so it's easy to "scrub the hub"
3. Corner grips for suction tubing and more
4. The adjustable support post lets you position it just where you need it
5. The anesthesia circuit clips keep your circuit fixed and firmly in place
6. Patented clamp attaches to OR bed, IV pole, or inpatient bed

**Anestand can help health systems save lives,  
lower costs and improve compliance.  
All at the same time.**





Scan to learn more

**The first and only anesthesia stand to safely and securely keep all supplies where they're needed can also be one of your best weapons in combatting infections.**

### **Anestand Facts**

- A small investment could help in compliance, save time and quickly pay for itself.
- By having supplies and medications organized and easily available with IV port ready, Anestand aids in patient safety.
- The bottom line: Anestand could help prevent costly infections and save money by helping anesthesia providers save OR time.

**Of all the many safety and security features Anestand offers, perhaps the single greatest is how it helps health systems to lower HAIs.**

- Having IV port ready for site alcohol scrub makes disinfection (and compliance) easier — Anestand keeps your IV port ready for scrubbing.
- By keeping medications immediately available and infusion ports clearly defined Anestand provides additional safety.
- By keeping supplies immediately available, Anestand helps make RSI (rapid sequence inductions) safe.
- During COVID, Anestand assists anesthesia providers as an extra pair of hands, while other healthcare personnel are being asked to keep distance from the airway.

**Anestand can provide value in safety, compliance, and convenience.**

**As soon as you use it, you'll wonder how you ever managed without it.**

- CDC, JCAHO, SHEA, and other organizations require cleaning the port with every injection. Anestand helps your facility and providers comply with these regulations.
- Patients with an HAI stay 19 additional days, or \$11.2 million in HAI treatment and/or lost revenue. For a small fixed cost your facility could get a large return on investment.
- Keeping supplies organized on a tray may prevent dropped items and reduce replacement costs. Anestand makes it far easier to keep supplies organized and right at hand.



# GUIDEBlade™

Precise, on the wire incisions for guidewire assisted catheter insertion

**GuideBlade—Precision Incisions, Every Time!**

## Evaluation of the efficacy and safety of a novel dermatotomy device for central venous cannulation

Worasak Keeyapaj, MD., Albert T. Cheung, MD.

Department of Anesthesiology, Perioperative and Pain Medicine, Stanford University, Stanford, California

Outcome	Number (%)
<b>Primary outcome:</b> Successful lines insertion without additional instrument	188 (100%)
Wire damage	0 (0%)
Operator injury	0 (0%)



## About the GuideBlade™

The GuideBlade™ is the only disposable scalpel on the market engineered specifically to streamline over the wire vascular access placement. The GuideBlade™ allows clinicians to perfectly align the dermatotomy incision made during vascular access placement with the in situ guidewire. This alignment prevents residual skin tags and eliminates the need for dermatotomy revisions.

## GuideBlade™ Clinical Benefits

- Less insertion site bleeding.
- Maintenance of insertion site skin integrity.
- Decreased risk of catheter bacterial colonization.
- Decreased risk of wire related complications (wire fracture and wire kinking).
- Decreased scarring/improved cosmetic result.

## GuideBlade™ In Practice

- Appropriate for use with guidewires up to and including 0.035" diameter.
- Useful for placement in a wide range of venous access devices including central venous catheters, PICC lines, venous introducer sheaths, hemodialysis catheters, and percutaneous ECMO cannulas.
- Particularly useful in procedures placed while wearing protective faceshields where condensation and glare can obscure visualization of the guidewire insertion site.
- Particularly useful in obese or bleeding patients where wire entry site can be obscured by excess soft tissue or entry site bleeding.
- Increased benefit when used by trainees and less experienced practitioners.

## The GuideBlade™ Value Proposition

- Decreases the need to replace kit components (like guidewires, dilators or catheters) that are easily damaged by residual skin tags often seen with standard disposable scalpel use.
- Decreases insertion site bleeding, thereby decreases the need for frequent dressing changes.
- Maintains skin integrity and thereby decreases risk of catheter colonization and associated infection.

Misaligned skin incisions make passage of dilators and catheters difficult and can lead to guidewire kinking and fracture. Revision of incisions may increase the risk of insertion site bleeding, scarring and catheter associated bloodstream infection.





Catalog 5.0