

B.Well

Caring for everyone



PRODUCT CATALOGUE

Home health
care solutions

 PRO  MED  TECHNO

BLOOD PRESSURE

RESPIRATORY

FEVER

ORAL HYGIENE

B.Well FAMILY

**A healthy family
is the basis of a happy life.
Nothing brings as much joy
as good health of one's family.**

B.Well Swiss sincerely shares people's desire to care for their loved ones. For this purpose B.Well Swiss has developed a wide range of medical devices meeting the needs of every member of the family. Helpless babies, active yet careless teenagers, elderly parents – they all need care and attention. B.Well Swiss provides diverse devices for all members of the family: blood pressure monitors, nebulizers, thermometers, irrigators, as well as series of devices for children under 3 years of age.

B. Well Swiss products guarantee perfect monitoring of the entire family's health, meaning protection and demonstrating one's care and concern.

B.Well Swiss advocates family values and caring for those in need of help, and bases its development on the health and well-being of all members of society.

B.Well product family intends to bring confidence and peace to every household, making more and more families happy!





High quality provider

CARING FOR EVERYONE



B.Well LINES

We made your choice easier by dividing all our products into 3 lines: MED, PRO and TECHNO. Now everyone can find the necessary product matching his or her needs.



PRO medical devices are easy to use, professional high-quality devices intended for prophylaxis, treatment and control of your family's health.



MED lines comprises devices with advanced functions. They enable you to apply up-to-date diagnostic methods, prophylaxis and treatment as well as conduct overall health monitoring for your family.



TECHNO devices have the advantages of state of the art technologies. That is why they are highly efficient devices that offer breakthrough opportunities enabling you to care for every member of your family.



Product was designed, verified and supervised by B.Well Swiss AG, Switzerland. Applied Swiss control system assures high quality, endurance and safety of all B.Well products.

SWISS HEALTH TECHNOLOGIES

The company was established in England in 2004 with the aim of producing high quality medical devices. In 2012, the head office was moved to Switzerland. This country symbolizes peacefulness and well-being, which are characteristics of a close-knit, healthy family. The true quality of B.Well Swiss's range of products is provided by working with highly skilled precision devices specialists, as well as cooperating with the best doctors and medical institutions.

MISSION

Creating products that help people take care and protect the health of their loved ones and making millions of families happy.

CARE

B.Well Swiss products help people take care of their families and change their lives for the better.

PROTECTION

The quality of B.Well Swiss products ensures control and protection of the family's health.

RELIABILITY

Thanks to excellent management, the production process of B.Well Swiss's products is as precise as Swiss watches.

EASY TO USE

B.Well Swiss products are easy to use, they don't require special skills or expertise, and provide high efficiency.

AVAILABILITY

Openness and focus on consumers and partners make it possible to create medical products available to everyone.

B.Well Swiss pays close attention to the quality control of its products. That is why the quality management system of the company is certified to comply with both European and international standards.



QUALITY CERTIFICATES



General Manager
of B.Well Swiss AG
Bernadette Maeder-Holenstein

The quality management system of B.Well Swiss AG is certified to comply with the international norm ISO 13485:2016. B.Well Swiss products comply with the Medical Device Directive (MDD) 93/42/EEC or Medical Device Regulation (EU) 2017/745.

The implementation of these standards confirms the status of B.Well as a reliable partner who offers a long-term and mutually beneficial relationship for costumers.



High quality provider

B.Well 

BLOOD PRESSURE MONITORS



BLOOD PRESSURE MONITORING

guarantees good health and wellbeing. Hypertension control will prevent dangerous health problems – heart attacks and strokes, cardiac and renal failures. Blood pressure measurement is simple and comfortable with B.Well's digital blood pressure monitors.

B.Well blood pressure monitors represent reliable protection for your family. The company's engineers have created blood pressure monitors for the needs of different people: The range of B.Well devices include both affordable blood pressure monitors and models with extended functionality.

Upper arm automatic blood pressure monitors are designed to provide comfortable and accurate measurements at home and at clinics. Wrist devices are comfortable to carry and easy to use on trips or at the office.

Recent developments allow us not only to correctly measure the pressure, but also to point out the signs of diseases such as hypertension and to identify early symptoms of arrhythmia.

B.Well blood pressure models are distinguished by their technological effectiveness, stylish design, and easy usage. B.Well blood pressure monitors are accurate, reliable, and available for everyone.

SWISS HEALTH TECHNOLOGIES



An advanced measurement accuracy is achieved by the automatic intellectual analysis of three successive measurements.



Measurement during inflation: individual level of air inflation for fast and gentle measurement.



Considering WHO's recommendation the blood pressure monitor calculates an average of the last 3 measurements.



Micro USB port for power supply: convenient and always at hand.



LCD screen color backlight warns against high blood pressure level.



Display backlight for comfortable use in any light environment.




Measurement during deflation: comfortable measurement, accurate result



Serves to detect pulse irregularities at an early stage



Blood pressure indication scale helps to interpret measurement values



Easy handling with only one button



Anatomic cuff guarantees comfortable and accurate measurement



RECOMMENDATION

Table for classifying blood-pressure values (unit: mmHg) according to European Society of Hypertension (ESH)

RANGE	SYSTOLIC BLOOD PRESSURE	DIASTOLIC BLOOD PRESSURE	RECOMMENDATION
Grade 3: severe hypertension	higher or equal to 180	higher or equal to 110	Urgently seek medical advice!
Grade 2: moderate hypertension	160-179	100-109	Consult your doctor immediately
Grade 1: mild hypertension	140-159	90-99	Consult your doctor
High intended	130-139	85-89	Consult your doctor
Normal	Lower than 130	Lower than 85	Self-check
Optimal	Lower than 120	Lower than 80	Self-check

B.Well 

MED-55 Automatic

Intelligent measurement

- Advanced blood pressure control: 60 memories for each of 2 users with date and time
- Anatomic cuff for an upper arm circumference of 22-42 cm

 MED


micro USB



Memory



Pulse Arrhythmia Detection



Comfortable anatomic cuff



Cuff check indicator: the device shows if the cuff is properly applied



Up-to-date algorithm of measurement

Intellect Active | Measurement during inflation

Individual level of air inflation for fast and gentle measurement



Micro USB port for power supply: convenient and always at hand

MED-55 blood pressure monitor provides an advanced measurement accuracy due to the 3check technology. Because blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement. You can choose between the 3check mode with a triple measurement and the single measurement mode.

Measurement 1
15-second pause
Measurement 2
15-second pause
Measurement 3
Analysis of 3 measurements
Result is as accurate as at a doctor's office



Unique technology

- Three consecutive measurements in the automatic mode
- Analysis of the measurement data based on a special algorithm
- Maximized accuracy of blood pressure calculation
- Where necessary, the device automatically performs a fourth, additional measurement



Colored display backlight (green, yellow or red) indicates the level of blood pressure



MED-53 Automatic

Gentle accurate measurement

- Intellect Active – fast and gentle measurement
- Micro USB port for power supply
- Large display with backlight
- Cuff check indicator: the device shows if the cuff is properly applied
- 90 memories with date and time
- Anatomic cuff for an upper arm circumference of 22-42 cm



Display
backlight



Memory



Pulse Arrhythmia
Detection



Comfortable
anatomic cuff

♥ MED



MED-57 Wrist

Gentle accurate measurement

- Intellect Active – fast and gentle measurement
- Average of the last 3 measurements
- Large display with backlight
- 30 memories with date and time
- For wrist circumference of 13,5-21,5 cm



Display
backlight



Memory



Pulse Arrhythmia
Detection



Average of
3 measurements

♥ MED

B.Well 


PRO-35 Automatic

High accuracy and measurement evaluation

- Intellect Classic – comfortable measuring, accurate result
- Memory of up to 30 measurements
- Anatomic cuff for an upper arm circumference of 22-42 cm
- Blood pressure indication scale


 PRO


PRO-33 Automatic

High-accuracy measurements

- Intellect Classic – comfortable measuring, accurate result
- Last reading memory
- Anatomic cuff for an upper arm circumference of 22-42 cm


 PRO

PRO-30 Semi-automatic

Affordable blood pressure control

- Easy handling with only one button
- Anatomic M - size cuff for an upper arm circumference of 22-32 cm
- Comfortable PVC bulb with quick air release valve

PRO



PRO-39 Wrist

High-accuracy measurements

- Compact and portable for travel, home and office
- Intellect Classic – comfortable measuring, accurate result
- Last reading memory
- For wrist circumference of 14-19.5 cm

PRO



B.Well 


MODEL	Size	Upper arm circumference
WA-C-ML	M-L	22-42 cm
WA-C-M	M	22-32 cm

ANATOMIC CUFFS

Comfortable and accurate measurements

- Fan-shaped for the highest level of comfort
- Available in M and M-L sizes
- Suitable for all common arm sizes

Anatomic cuff provides comfortable and accurate measurements. Conical form is convenient to wear and cuff wraps easily around the arm.

M-L size cuffs are suitable for all, thus allowing them to be used by all family members.

AD-53/55 Micro USB mains adapter

(for models MED-53, MED-55)



- 100-240 V, 50/60 Hz AC
- 6 V, 600 mA DC

ASSA51E mains adapter

(for models PRO-33, PRO-35, PRO-36)



- 100-240 V, 50/60Hz AC
- 6 V, 600 mA DC



RECOMMENDATION

For taking reliable measurements

- Cuff is as important as the blood pressure monitor
- The conical shape provides a better fitting and accuracy
- Using the wrong size cuff is a major cause for erroneous blood pressure measurement

ANEROID BLOOD PRESSURE KITS AND STETHOSCOPES



SPHYGMOMANOMETERS & STETHOSCOPES

are included in B.Well's product range for specialists.

B.Well aneroid sphygmomanometers have gained trust and popularity among healthcare specialists in many countries. These devices meet high accuracy and reliability standards. B.Well's range of aneroid devices provide a choice of home use devices for everyone who possesses the necessary skills and who prefers

the conventional measurement method. B.Well sphygmomanometers are accurate, convenient, and durable.

The main advantage of B.Well stethoscopes is their excellent sound transmission. At the same time they are affordable, yet of high quality and are available in different colors, thus giving you a range of suitable options.



High quality provider

B.Well 

Bulb with a dust filter
for accurate
measurement

MED-62

Stethoscope included

Reliable and durable

- Cuff with a metal ring for professional and home use
- Soft storage bag



Universal cuff
with a ring



Long-lasting
performance



Bulb
dust filter



Durable
bladder



Universal cuff for
different upper arm
size



Long-lasting
metal gauge



Printing:
steady
against
rubbing



Universal size 22-42 cm, suits most
adults



Rounded shape: easy to put on





MED-63

Built-in stethoscope

Reliable and durable

- Cuff with metal ring for professional and home use
- Built-in stethoscope for comfortable usage
- Storage bag



♥ MED



MED-61

Professional

Essential to have

- Cuff without ring for professional use
- Storage bag



♥ MED

B.Well 


PRO-60

Stethoscope included

Reliable and durable

- Cuff with a metal ring for professional and home use
- Stethoscope included – ready for use


 PRO


WS-1

Single head

Lightweight

- Soft cone-shaped ear tips
- Lightweighted and compact
- Comfortable chrome-plated binaural
- Durable one-piece molded PVC Y-tube
- Plastic storage bag


 bordeaux

 blue


 MED



WS-2 Dual head

Versatile and affordable

- Bell with a non-chilling ring provides comfort
- Soft cone-shaped ear tips
- Comfortable chrome-plated binaural
- Durable one-piece molded PVC Y-tube
- Plastic storage bag



red



lilac

 MED


WS-3 Sprague Rappaport Type Stethoscope

Excellent acoustics

- Multifunctionality in one unit
- Double tubes provide excellent acoustics
- Soft cone-shaped ear tips
- Furnished with 3 different sized bell chest pieces
- 3 different spare ear tip sizes
- Spare plastic membranes in 2 different sizes
- Plastic storage bag



red



lilac



blue



dark grey

 MED

PRODUCTS FOR RESPIRATORY THERAPY



INHALATION THERAPY

is the most effective approach for treating and preventing respiratory diseases. It delivers the medication directly to the target organ and ensures rapid recovery without any adverse effects on other organs.

B.Well medical inhalers and nebulizers have proved to be reliable and effective. At the same time the company is working to improve its products using up-to-date technological achievements. Design and modification processes involve competent engineers and

pulmonologists, which makes it possible to achieve fruitful and simple solutions to maintain the health of loved ones.

B.Well inhalers come in many types for treating the upper and lower respiratory tract and the range includes a compressor, mesh, and steam nebulizer, thus letting everyone choose their treatment approach.

B.Well nebulizers are designed for effective treatment and convenience for every member of the family!

NEW!

MED-120

Compact compressor nebulizer

Effective and convenient

- Reliable compressor with long-lasting performance
- All the drugs can be used for inhalation therapy
- Saves drugs – low residual volume of medicine
- Nebulizer holder for convenient use
- Soft bag included for convenient transportation and storage.


 MED

RESPIRATORY

BLOOD PRESSURE



Micro USB



Micro USB port for power supply: convenient and always at hand



Quiet operation ($\leq 45\text{dB}$) and ultra compact size – easy to use, carry and store



Two masks – comfortable for children and adults

B.Well MED-120 compressor nebulizer ensures effective therapy in the lower respiratory tracts. Particle size ($2,9\text{ }\mu\text{m}$) and high respirable fraction (71%) allow medication to penetrate deep into the lungs for fast healing.

Micro USB is easy and convenient for use in any daily environment
Compressor nebulizer MED-120 conforms to European standard for nebulizers EN13544-1.



High quality provider


B.Well 



TH-134 Mesh nebulizer

Effective and innovative

- Average particle size < 5 μm – effective therapy of lower respiratory tract
- Portable (137 g) – easy carriage and prophylactic treatment of asthma
- Noiseless while operating at an angle of up to 45° – possible to perform inhalation while sleeping or in comfortable position




TECHNO



TH-135 Mesh nebulizer

Special child version

- Two masks – comfortable for children and adults
- AC adapter included
- Operates on batteries and mains power supplies – easy usage at home and on trips




TECHNO



MED-121 Compressor nebulizer

Comfortable and effective

- Quiet and lightweight
- Innovative nebulizer «Family» (MMAD $\leq 3,0$ micron) – effective lower respiratory tract therapy
- Nebulization rate selection for children and adults
- Two masks – comfortable for children and adults




♥ MED



MED-125 Compressor nebulizer

Special child version

- Ideal for children starting from birth – including child and infant masks



♥ MED


B.Well 


PRO-110

Compressor nebulizer

Professional and fast treatment

- Modern nebulizer «Basic» (MMAD ~ 3,16 micron) – effective lower respiratory tract therapy
- Nebulization speed 0,4 ml/min – intense and rapid inhalation
- Two masks – comfortable for children and adults
- Reliable compressor for long-lasting performance
- Nose piece helps to cure stuffy nose


 PRO


PRO-115

Compressor nebulizer

Effective and safe treatment

- Toy nebulizer – playful therapy for children!
- Three masks – comfort for the whole family
- Modern nebulizer «Basic» (MMAD ~ 3,16 micron) – effective lower respiratory tract therapy
- Reliable compressor for long-lasting performance
- Nose piece helps to cure stuffy nose


 PRO





PRO-118 Steam inhaler

Effective cold prevention

- Two temperature modes – comfortable for children and adults
- Cosmetic treatment mask – gentle care for skin

PRO



Effective prevention
against cold,
pharyngitis,
and laryngitis



Possible to use decoctions,
herbal extracts, essential oils,
mineral and salt water
for treatment and care



Mild steam temperature
of 43°C benefits and does not
burn skin

Steam inhaler generates steam with particles bigger than 10 microns, which effectively treat the upper respiratory track.

Thanks to this device, traditional inhalation becomes more comfortable and easier and safer.



RECOMMENDATION

Water solution preparation:

1 table spoon of crushed dry herb per 200 ml hot water. Prepare and thoroughly filter (prepared solution should be clear and should not contain visible herbal particles).


Spirit solution preparation:

20 drops of spirit solution (a half a teaspoon) per 100 ml saline solution.



High quality provider

TECHNICAL DATA

LINE	 PRO	 PRO	 PRO	 MED	 MED	 MED	 TECHNO
MODEL	PRO-118	PRO-110	PRO-115	MED-120	MED-121	MED-125	TH-134 & 135
							
TYPE	steam inhaler	compressor		compressor	compressor		mesh
SIZE (mm)	235 x 90 x 260	137 x 173 x 96	224 x 112 x 170	100 x 54 x 51	102 x 140 x 61		45 x 54 x 122
WEIGHT (without accessories and batteries)	~ 0,560 kg	1,345 kg	1,51 kg	0,235 kg	0,303 kg		0,137 kg
MEDICATION CAPACITY	80 ml	2-8 ml		2-6 ml	2-5 ml		2-8 ml
PARTICLE SIZE (microns μm)	> 10 μm	MMAD ~ 3,16 μm % of particle size < 5 μm (FDP) $\geq 70\%$		MMAD $\leq 2.9 \mu\text{m}$ % of particle size < 5 μm (FDP) $\geq 71\%$	MMAD $\leq 3.0 \mu\text{m}$ % of particle size < 5 μm (FDP) $\geq 70\%$		< 5 μm
NEBULIZATION RATE	~ 10,6 ml/min	0,4 ml/min		0,3 ml/min	Fully open valve: ≥ 0.30 ml/min Fully closed valve: ≥ 0.08 ml/min		0,2 – 1,0 ml/min

COMPLETE ACCESSORIES SET

MODEL / ACCESSORIES SET	PRO-118	PRO-110	PRO-115	MED-120	MED-121	MED-125	TH-134	TH-135
CHILD MASK	-	✓	✓	✓	✓	✓	-	✓
INFANT MASK	-	-	✓	-	-	✓	-	-
ADULT MASK	-	✓	✓	✓	✓	-	-	✓
MOUTHPIECE	-	✓	✓	✓	✓	✓	✓	✓
NOSEPIECE	-	✓	✓	-	-	-	-	-
INHALATION MASK	✓	-	-	-	-	-	-	-
COSMETIC TREATMENT MASK	✓	-	-	-	-	-	-	-
AIR TUBE	-	✓	✓	✓	✓	✓	-	-
AIR FILTERS (5 pcs)	-	✓	✓	✓	✓	✓	-	-
MAINS SUPPLY	✓	✓	✓	✓	✓	✓	-	-
BATTERIES (2 pcs)	-	-	-	-	-	-	✓	✓
AC ADAPTER	-	-	-	✓	✓	✓	-	✓
STORAGE BAG	-	✓	✓	✓	✓	✓	✓	✓
STICKERS	-	-	✓	-	-	✓	-	-

INHALER CHOOSING ALGORITHM



COMMON MEDICINE AND INFUSIONS FOR RESPIRATORY THERAPY

MEDICINE TYPE	NEBULIZER		INHALER
	COMPRESSOR	MESH	STEAM INHALER
Corticosteroids for inhalation	+	+	—
Mucolytics	+	+	—
Antibiotics, antiseptics	+	+	—
Mast cell stabilizer, antihistamine	+	+	—
Anticholinergics	+	+	—
beta2-adrenoceptor agonists	+	+	—
Combined (Composite) medicine	+	+	—
Normal saline	+	+	+
Essential oil	—	—	+
Herbal infusions	—	—	+

DIGITAL AND INFRARED THERMOMETERS



THERMOMETER

is an essential device for any family's home. High accuracy, safety, and convenience are the advantages of B.Well Swiss infrared and digital thermometers. The company offers devices which are infant-friendly, handy for parents, and comfortable for the whole family.

Digital and infrared thermometers are safe to handle and do not contain glass or mercury, as opposed to mercury thermometers.


B.Well thermometers comply with international quality and safety standards.

The main advantage of infrared thermometers is a fast measurement within a second. It measures the body temperature without touching the skin and therefore does not disturb the infant. Infrared thermometers can be used to take body, bath water, room air, etc. temperature.




Multifunctional


B.Well infrared thermometers can be used to take body, bath water, room air, etc. temperature. They are universal devices which can be used instead of other thermometers.




B.Well infrared thermometers incorporate a wide measurement range feature




Automatic storage of the last 10 or 32 measurements – enabling efficient tracking of past measurements



Backlit display enables one to see temperature readings even at night. Illumination has a soft luminescence, which guarantees comfort for the eyes and great visibility



Flexible tip adjusts to body form for a comfortable measuring



Thermometer is completely waterproof

NORMAL TEMPERATURE VALUES

MEASURING METHOD	NORMAL RANGE
Axillary	35.2-36.8 °C
Oral	35.7-37.3 °C
Rectal	36.2-37.7 °C
Ear	35.5-37.5 °C
Forehead	35.5-37.3 °C



RECOMMENDATION

It is possible to take axillary, oral, and rectal measurements.

To receive reliable results, we recommend taking oral or rectal measurements for infants and small children. These are very reliable because the thermometer sensor does not measure the air temperature.



High quality provider

B.Well 

NEW!





WF-4000

Professional measurement

- Versatile 3 in 1: Measures the temperature of a human body, water, air and objects
- A special signal will inform that temperature is elevated
- Made of safe materials: thermometer is made of safe plastic and contains no mercury or glass

TECHNO



Contact-free measurement: exact result in only 1 second



Comfortable measurement at night: a large backlit display, sound signal can be turned off



Memory for 32 measurements allows to follow temperature dynamics



Ergonomic design ensures comfort and ease of use

Non-contact infrared thermometer WF-4000 is swichable from °C to °F scales.



WF-4000 measures not only the temperature of the human body, but also the tempareture of water, air and objects.

WF-5000

Contact-free and accurate measurements

- Measures temperature of objects, air, and liquids
- Provides high measurement accuracy due to indication of correct distance
- Quick Action – measurements within one second
- User-friendly – easy to use thanks to display icons
- Soft bag for convenient storage

TECHNO



Distance sensor indicates the right position for an accurate measurement



Illuminated display for easy measuring at night



10 data memory function tracks past temperature measurements



Fever alarm with special beep

Thermometer saves parents's time by measuring kids temperature within one second. WF-5000 measures temperature without touching and

therefore does not disturb a sleeping baby. Backlit display will allow to see temperature measurements at night. WF-5000 is swichable from °C to °F scales.





High quality provider


B.Well 


WF-1000

Instant and accurate measurements

- Measurements possible both in the ear canal and on the forehead
- Measures temperature of objects, air, and liquids
- Accurate measurement within seconds





Waterproof sensor is reliable and easy to clean



User friendly: Icons on the display make measuring clear and easy



Instant measurement both in the ear canal and on the temple

It takes seconds to measure and display the temperature. Ideal for children and infants –

quick temperature measurement without disturbing the sleeping child or infant.



RECOMMENDATION


- Before the measurement, patients as well as the thermometer should stay in a steady state condition for at least 30 minutes
- Allow one minute between successive measurements as slight variations may occur if measurements are taken over a short period of time
- Be sure to clean the thermometer lens each time after usage

WT-06 flex

Accurate and fast

- Flexible tip for gentle touch
- Beeper and Fever Alarm – acoustic signal if temperature is high
- Case for convenient storage

PRO



Children will enjoy the funny duck or rabbit



Made of safe materials



Shows the last measurement

WT-06 flex thermometer has a flexible tip and provides a gentle skin touch. Temperature measuring time is about 10 seconds, which does not distract the infant for long.

Due to its kids design and bright colors, child will enjoy this thermometer. The bright duck or cheerful rabbit will become the kid's best friend and turn any measurement into a game.



RECOMMENDATION

- For receiving exact figures it is necessary for the tip to have a solid contact with the body. If the tip is not in close contact with the skin the unit beeps earlier and the measured temperature is inaccurate. Keep mouth closed when taking an oral temperature. When taking axillary the temperature keep the arm close to the body so the probe tip is well covered by the arm.



High quality provider

B.Well 



WT-04 standard

Accurate and comfortable

- Flexible tip for gentle touch
- Anti-skid panels – comfortable to hold
- Displays most recent measurement
- Case for convenient storage




PRO



WT-03 base

Accurate and reliable

- Displays most recent measurement
- Turns off automatically for energy efficiency
- Case for convenient storage



PRO





WT-09 quick


Comfort for babies

- Pacifier thermometer will calm your baby
- Displays most recent measurement
- Made of safe materials
- Case for convenient storage

♥ MED










CHOOSING A THERMOMETER



High quality provider

TECHNICAL DATA

			DESCRIPTION	MEASUREMENT RANGE	TIP	MEASUREMENT TIME	WATERPROOF	MEASURING OPPORTUNITIES	CONTACT/ NON-CONTACT
PRO	WT-03 base		For the whole family	32-43°C	not flexible	60 sec	yes	Body temperature	Contact
	WT-04 standard		Comfortable model with a flexible tip	32-43°C	flexible	30 sec	yes	Body temperature	Contact
	WT-06 flex		Bright and fast: the best choice for a kid	32-43°C	flexible	6-10 sec	yes	Body temperature	Contact
MED	WT-09 quick		Pacifier thermometer for infants	32-43°C	pacifier	1,5-3 min	100% waterproof	Body temperature	Contact
	WF-1000		Infrared ear/forehead thermometer for a reasonable price	10-50°C		from 1 sec	-	Body, object, and ambient temperature	Contact
TECHNO	WF-4000		Advanced non-contact infrared thermometer	0-110°C (32 °F - 230 °F)		from 1 sec	-	Body, object, and ambient temperature Switchable °C to °F	Non-contact
TECHNO	WF-5000		Advanced non-contact infrared thermometer	0-100°C (32 °F - 212 °F)		from 1 sec	-	Body, object, and ambient temperature Switchable °C to °F	Non-contact



WC-150

Nasal aspirator «CleanNose»

PRO

Clean nose, easy breathing

- Cleans baby's nose secretions in a few seconds, prevents infections, improves the baby's health, and speeds up their recovery
- Compressor creates ideal pressure for effective and gentle nasal cleaning
- Easy handling with only one button
- Soft storage bag



Gentle and safe



Hygienic



Easy handling



For children from birth



Soft silicone tips prevent any injuries of the baby's soft skin



Special design ensures easy accessory detachment and cleaning



12 well-known children's tunes are stored on the memory, which will distract the baby while operating

Tips and purpose

The large tip is recommended for use when the mucus is in a more solid form



The narrow tip is recommended for use when the mucus is loose and in a liquid form



| ORAL HYGIENE

B.Well oral irrigators are effective devices for daily oral hygiene to be used by every member of the family to maintain dental health and beauty!

There are hard-to-reach spots in the oral cavity in which microbes accumulate. A tooth brush may remove only half of bacteria. Incomplete cleaning of the teeth results in caries and gingival diseases.

That is why dentists constantly remind people that an irrigator is to be used daily in addition to a tooth brush!




NEW!

WI-912 Oral irrigator

Modern care for teeth
and gums

- Maximum power of 720 kPa for complete hygiene
- Nozzles stand
- 360° nozzle rotation
- Wall mounting kit

TECHNO



Convenient bag
for storage and
transportation



Non-contact induction
charger combined
with irrigator and nozzles
stand



Wall mounting
kit - saves place
in the bathroom



Why should you use an irrigator:

- Prevention of caries and gingival diseases
- Removal of plaque which stimulates bacteria formation
- Prevention of dental tartar
- Prevention of bad breath
- Appropriate care if you have dental implants, crowns, bridges, braces, and retainers
- Complete and effective interdental cleaning
- Therapeutic gingival massage



High quality provider

NEW!

WI-922

Oral irrigator

Professional care of your teeth and gums

- Pulse mode for maximum efficiency
- Ten modes for smooth flow adjustment and for individual intensity selection
- 360° nozzle rotation for complete hygiene
- Convenient holder for the device handle
- Optimal water tank capacity 600 ml
- Anti-slip stable feet
- Convenient cover case for nozzles storage



7
NOZZLES
in a set



360° nozzle rotation



Family








Water tank capacity 600 ml



Adjustable power

Nozzles types and purpose:

Standard nozzle for oral care	Nozzle for tongue cleaning and fresh breath	Orthodontic nozzle for braces	Periodontal nozzle for gums	Nozzle for implants, dental bridges, crowns and veneers
				



WI-911 Portable oral irrigator

A beautiful smile and
fresh breath

PRO

- AC adapter included in the set – device can be charged anywhere
- Easy water refilling
- Low noise level
- Two replacement nozzles in different colors – ideal for families or couples
- Auto shut-off after 2 minutes – device does not discharge if accidentally turned on



3 modes
of operation



Optimal power




Compact



Pouch as a gift



Cordless. Equipped with
modern Li-ion battery –
possible to use for one
month on one charge



Optimal power
of 620 kPa for
effective hygiene



Three modes of operation
for comfortable usage:
Normal – effective hygiene
Soft – for sensitive teeth
Pulse – therapeutic gingival massage

Set content:

Portable
oral irrigator



Nozzles
(2 pcs)



Water
tank
(150 ml)



AC adapter



High quality provider





B.Well Swiss AG

Bahnhofstrasse 24, 9443 Widnau,
Switzerland

info@bwell-swiss.ch
www.bwell-swiss.ch

